



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 03/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101306

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1687-01	HR-06-031119-A	1	1.13	9.55	8.28	84.9
K1687-03	HR-06-031119-B	2	1.12	9.94	8.47	83.3
K1687-05	HR-06-031119-C	3	1.15	9.74	8.39	84.3
K1693-05	OR-02-031219-A	96	1.15	9.95	9.22	91.7
K1693-07	OR-02-031219-B	97	1.12	9.65	8.95	91.8
K1697-01	SU-03-031119-A	4	1.16	9.74	8.42	84.6
K1697-03	SU-03-031119-B	5	1.18	9.72	8.78	89.0
K1697-05	SU-04-031119-A	6	1.18	9.97	9.37	93.2
K1697-07	SU-04-031119-B	7	1.18	9.74	8.89	90.1
K1807-02	TP-1	8	1.17	9.98	8.97	88.5
K1807-04	TP-2	9	1.16	9.8	8.42	84.0
K1807-06	TP-3	10	1.15	9.58	8.24	84.1
K1807-08	TP-4	11	1.12	9.62	8.69	89.1
K1807-10	TP-5	12	1.14	9.77	8.63	86.8
K1807-12	TP-6	13	1.19	9.68	8.58	87.0
K1807-14	TP-7	14	1.15	9.95	8.94	88.5
K1807-16	TP-8	15	1.18	9.94	8.8	87.0
K1807-18	TP-9	16	1.19	9.86	8.63	85.8
K1807-20	TP-10	17	1.16	9.77	9.00	91.1
K1807-22	TP-18	18	1.15	9.83	9.02	90.7
K1807-24	TP-17	19	1.17	9.81	8.52	85.1
K1807-26	TP-16	20	1.13	9.77	8.67	87.3
K1807-28	TP-15	21	1.1	9.81	8.79	88.3
K1807-30	TP-14	22	1.17	9.59	8.53	87.4
K1807-32	TP-13	23	1.15	9.51	8.49	87.8
K1807-34	TP-12	24	1.16	9.55	8.91	92.4
K1807-36	TP-11	25	1.15	9.94	8.78	86.8
K1807-38	TP-19	26	1.18	9.89	8.8	87.5
K1807-40	TP-20	27	1.16	9.85	8.61	85.7
K1807-42	TP-21	28	1.18	9.98	9.28	92.0
K1807-44	TP-22	29	1.11	9.75	8.72	88.1
K1807-46	TP-23	30	1.12	9.89	8.36	82.6
K1824-02	TP-24	31	1.13	9.55	7.89	80.3
K1824-04	TP-25	32	1.15	9.65	8.38	85.1
K1824-06	TP-26	33	1.15	9.95	8.61	84.8
K1824-08	TP-27	34	1.16	9.54	8.44	86.9
K1824-10	TP-28	35	1.15	9.88	8.86	88.3
K1824-12	TP-29	36	1.1	9.66	8.69	88.7
K1824-14	P-1	37	1.14	9.64	8.71	89.1



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 03/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101306

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1824-16	P-2	38	1.13	9.71	8.49	85.8
K1824-18	P-3	39	1.19	9.82	8.59	85.7
K1824-20	P-4	40	1.14	9.79	8.44	84.4
K1826-01	B-1	41	1.18	9.64	8.26	83.7
K1826-02	B-2	42	1.16	9.84	8.63	86.1
K1826-03	B-3	43	1.16	9.95	8.3	81.2
K1830-01	BALLFIELD-SS-03072019	44	1.15	9.52	8.81	91.5
K1840-07	BP-VPB201-GW-640-642	45	1.13	9.96	2.05	10.4
K1843-01	GATE-SOIL-BOX	46	1.13	9.62	8.63	88.3
K1844-01	DRUM-2-4	47	1.00	2.00	2.00	100.0
K1844-09	TOTE-1-3	48	1.00	2.00	2.00	100.0
K1845-01	ETGI-294	49	1.00	2.00	2.00	100.0
K1846-01	TP-2-1	50	1.13	9.55	8.73	90.3
K1846-02	TP-2-2	51	1.16	9.75	8.72	88.0
K1846-03	TP-2	52	1.16	9.96	9.12	90.5
K1846-05	TP-1-S1	53	1.1	9.63	8.88	91.2
K1846-06	TP-1-S2	54	1.16	9.83	8.86	88.8
K1846-07	TP-1-S	55	1.16	9.83	9.01	90.5
K1846-09	TP-1-D1	56	1.17	9.96	8.76	86.3
K1846-10	TP-1-D2	57	1.16	9.98	9.05	89.5
K1846-11	TP-1-D	58	1.18	9.73	8.41	84.6
K1847-01	SUP-7B-G-1-031119	59	1.1	9.96	9.19	91.3
K1847-02	SUP-7B-G-2-031119	60	1.13	9.94	9.17	91.3
K1847-03	SUP-7B-G-3-031119	61	1.12	9.58	8.74	90.1
K1847-04	SUP-7B-G-4-031119	62	1.17	9.98	8.91	87.9
K1847-05	SUP-7B-G-5-031119	63	1.13	9.96	8.65	85.2
K1847-06	SUP-7B-G-6-031119	64	1.14	9.61	8.68	89.0
K1847-07	SUP-7B-G-7-031119	65	1.11	9.7	8.9	90.7
K1847-08	SUP-7B-G-8-031119	66	1.13	9.83	8.74	87.5
K1847-09	SUP-7B-G-9-031119	67	1.17	9.87	8.98	89.8
K1847-10	SUP-7B-G-10-031119	68	1.13	9.95	8.75	86.4
K1847-11	SUP-7B-G-11-031119	69	1.14	9.85	8.87	88.7
K1847-12	SUP-7B-G-12-031119	70	1.13	9.86	8.4	83.3
K1847-13	SUP-7B-G-13-031119	71	1.15	9.98	8.63	84.7
K1847-14	SUP-7B-C-1-031119	72	1.14	9.82	8.75	87.7
K1847-15	SUP-7B-C-2-031119	73	1.14	9.61	8.66	88.8
K1847-16	SUP-7B-C-3-031119	74	1.1	9.66	8.46	86.0
K1847-17	SUP-7B-C-4-031119	75	1.12	9.74	8.6	86.8
K1847-18	SUP-7B-C-5-031119	76	1.16	9.98	9.09	89.9



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 03/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101306

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1849-01	C1355-PC	77	1.00	2.00	2.00	100.0
K1849-03	H1308-PC	78	1.00	2.00	2.00	100.0
K1849-05	PUMP-PC	79	1.00	2.00	2.00	100.0
K1849-07	C135-CON	80	1.00	2.00	2.00	100.0
K1849-09	H1308-CON	81	1.00	2.00	2.00	100.0
K1849-11	PUMP-CON	82	1.00	2.00	2.00	100.0
K1851-01	390307016.1	83	1.13	9.52	1.64	6.1
K1851-02	390307011.2	84	1.16	9.98	1.41	2.8
K1851-03	390305716.0	85	1.13	9.92	4.66	40.2
K1851-04	390306745.0	86	1.15	9.9	2.81	19.0
K1851-05	390307026.0	87	1.16	9.86	1.55	4.5
K1857-01	390307113.0	88	1.16	9.96	5.52	49.5
K1857-02	390307733.1	89	1.13	9.98	2.54	15.9
K1857-03	390308025.0	90	1.13	9.96	1.57	5.0
K1861-01	BF-06-031119	91	1.12	9.57	8.77	90.5
K1861-02	BF-07-031119	92	1.16	9.97	9.18	91.0
K1861-03	BF-08-031119	93	1.2	9.64	8.74	89.3
K1861-04	BF-09-031119	94	1.14	9.71	8.97	91.4
K1861-05	BF-10-031119	95	1.15	9.59	9.00	93.0
K1863-01	MCH-CAULK	98	1.00	2.00	2.00	100.0
K1864-01	FLOOR	99	1.00	2.00	2.00	100.0
K1864-02	BRICK	100	1.00	2.00	2.00	100.0
K1864-03	BLOCK	101	1.00	2.00	2.00	100.0
K1865-01	DRUM-1-2	102	1.00	2.00	2.00	100.0
K1866-01	MR-MH-04-COMP	103	1.15	9.96	8.13	79.2
K1866-03	MR-MH-03-COMP	104	1.14	9.98	8.04	78.1
K1866-05	MR-MH-02-COMP	105	1.12	9.93	8.45	83.2
K1866-07	MR-MH-01-COMP	106	1.13	9.73	8.08	80.8
K1866-09	MR-MH-15-COMP	107	1.18	9.97	9.28	92.2
K1866-11	MR-MH-14-COMP	108	1.12	9.91	7.92	77.4
K1866-13	MR-MH-13-COMP	109	1.13	9.52	6.83	67.9
K1866-15	MR-MH-12-COMP	110	1.13	9.88	7.04	67.5
K1869-01	SP-1	111	1.14	9.92	8.81	87.4
K1869-02	SP-1-A	112	1.14	9.92	8.81	87.4
K1869-03	SP-1-B	113	1.14	9.92	8.81	87.4
K1869-04	SP-1-C	114	1.14	9.92	8.81	87.4
K1869-05	SP-1-D	115	1.14	9.92	8.81	87.4
K1869-06	SP-1-E	116	1.14	9.92	8.81	87.4
K1869-07	SP-1-F	117	1.14	9.92	8.81	87.4



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/13/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 03/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101306

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1869-08	SP-1	118	1.14	9.92	8.81	87.4
K1869-09	SP-2	119	1.11	9.86	8.73	87.1
K1869-10	SP-2-A	120	1.11	9.86	8.73	87.1
K1869-11	SP-2-B	121	1.11	9.86	8.73	87.1
K1869-12	SP-2-C	122	1.11	9.86	8.73	87.1
K1869-13	SP-2-D	123	1.11	9.86	8.73	87.1
K1869-14	SP-2-E	124	1.11	9.86	8.73	87.1
K1869-15	SP-2-F	125	1.11	9.86	8.73	87.1
K1869-16	SP-2	126	1.11	9.86	8.73	87.1
K1870-01	SP-1	127	1.14	9.94	8.79	86.9
K1870-02	SP-1-A	128	1.14	9.94	8.79	86.9
K1870-03	SP-1-B	129	1.14	9.94	8.79	86.9
K1870-05	SP-2	130	1.18	9.91	8.41	82.8
K1870-06	SP-2-A	131	1.18	9.91	8.41	82.8
K1870-07	SP-2-B	132	1.18	9.91	8.41	82.8
K1870-09	SP-3	133	1.14	9.68	8.51	86.3
K1870-10	SP-3-A	134	1.14	9.68	8.51	86.3
K1870-11	SP-3-B	135	1.14	9.68	8.51	86.3
K1870-13	ROW-CONCRETE-1	136	1.00	2.00	2.00	100.0
K1870-15	ROW-CONCRETE	137	1.00	2.00	2.00	100.0
K1871-01	TP-1	138	1.13	9.68	8.4	85.0
K1871-02	TP-1-A	139	1.13	9.68	8.4	85.0
K1871-03	TP-1-B	140	1.13	9.68	8.4	85.0
K1871-04	TP-1	141	1.13	9.68	8.4	85.0
K1871-05	TP-2	142	1.16	9.96	8.78	86.6
K1871-06	TP-2-A	143	1.16	9.96	8.78	86.6
K1871-07	TP-2-B	144	1.16	9.96	8.78	86.6
K1871-08	TP-2	145	1.16	9.96	8.78	86.6
K1873-01	371	146	1.1	9.82	8.46	84.4
K1873-02	371-A	147	1.1	9.82	8.46	84.4
K1873-03	371-B	148	1.1	9.82	8.46	84.4
K1873-04	371	149	1.1	9.82	8.46	84.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$