



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/14/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:30
In Date: 03/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 03/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101333

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1794-01	BF5N3	1	1.13	9.57	6.16	59.6
K1794-02	BF5N4	2	1.13	9.88	8.47	83.9
K1794-03	BF5N4MS	3	1.13	9.88	8.47	83.9
K1794-04	BF5N4MSD	4	1.13	9.88	8.47	83.9
K1794-05	BF5N5	5	1.16	9.98	6.06	55.6
K1794-06	BF5N6	6	1.14	9.98	5.76	52.3
K1794-07	BF5N7	7	1.15	9.96	7.17	68.3
K1794-08	BF5N8	8	1.17	9.92	6.5	60.9
K1794-09	BF5N9	9	1.16	9.97	6.34	58.8
K1794-10	BF5P0	10	1.13	9.57	5.94	57.0
K1794-11	VHBLK01	11	1.00	2.00	2.00	100.0
K1794-12	BF5N1	12	1.11	9.69	7.82	78.2
K1794-13	BF5N2	13	1.15	9.51	8.16	83.9
K1794-14	BF5P1	14	1.16	9.87	6.61	62.6
K1794-15	BF5P2	15	1.17	9.9	6.73	63.7
K1794-16	BF5P3	16	1.18	9.93	5.69	51.5
K1794-17	BF5P4	17	1.14	9.98	5.98	54.8
K1794-18	BF5P5	18	1.13	9.91	6.5	61.2
K1794-19	BF5P6	19	1.16	9.86	5.24	46.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$