



PERCENT SOLID

Supervisor: AMANDEEP
Analyst: jignesh
Date: 3/16/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:25
In Date: 03/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 03/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101393

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1936-01	BF0T8	1	1.14	9.82	6.94	66.8
K1936-02	BF0T9	2	1.12	9.93	6.98	66.5
K1936-03	BF0W0	3	1.15	9.74	7.07	68.9
K1936-04	BF0W1	4	1.13	9.71	7.78	77.5
K1936-05	BF0W2	5	1.16	9.95	7.4	71.0
K1936-06	BF0W3	6	1.12	9.69	7.81	78.1
K1936-07	BF0W4	7	1.16	9.64	7.14	70.5
K1936-08	BF0W5	8	1.1	9.96	7.72	74.7
K1936-09	BF0W6	9	1.18	9.72	7.75	76.9
K1936-10	BF0W8	10	1.13	9.97	7.49	71.9
K1936-11	BF0W9	11	1.16	9.64	6.97	68.5
K1936-12	BF0X0	12	1.13	9.82	7.00	67.5
K1936-13	BF0X1	13	1.15	9.86	7.11	68.4
K1936-14	BF0X2	14	1.18	9.98	7.59	72.8
K1936-15	BF0X3	15	1.15	9.91	7.4	71.3
K1936-16	BF0X4	16	1.14	9.94	7.97	77.6
K1936-17	BF0X5	17	1.19	9.63	8.46	86.1
K1936-18	BF0X6	18	1.13	9.71	8.59	86.9
K1936-19	VHBLK01	19	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$