



PERCENT SOLID

Supervisor: AMANDEEP
Analyst: jignesh
Date: 3/16/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:50
In Date: 03/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 03/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101394

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1939-02	BF1J4	1	1.19	9.98	7.64	73.4
K1939-04	VHBLK02	2	1.00	2.00	2.00	100.0
K1939-05	BF0S8	3	1.14	9.71	8.51	86.0
K1939-06	BF0S9	4	1.14	9.76	7.02	68.2
K1939-07	BF0T0	5	1.14	9.85	7.18	69.3
K1939-08	BF0T1	6	1.19	9.97	7.41	70.8
K1939-09	BF0T2	7	1.13	9.94	7.19	68.8
K1939-10	BF0T3	8	1.15	9.95	7.12	67.8
K1939-11	BF0T4	9	1.11	9.51	8.11	83.3
K1939-12	BF0T5	10	1.14	9.96	7.87	76.3
K1939-13	BF0T6	11	1.19	9.74	7.95	79.1
K1939-14	BF0T7	12	1.13	9.77	8.16	81.4
K1939-15	BF0X8	13	1.16	9.81	7.01	67.6
K1939-16	BF0X9	14	1.17	9.57	6.77	66.7
K1939-17	BF0Y0	15	1.15	9.5	6.69	66.3
K1939-18	BF0Y1	16	1.17	9.83	7.71	75.5
K1939-19	BF0Y2	17	1.14	9.79	8.07	80.1
K1939-20	BF0Y3	18	1.1	9.94	7.89	76.8
K1939-21	BF0Y4	19	1.16	9.96	7.98	77.5
K1939-22	BF0Y5	20	1.17	9.98	8.22	80.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$