



PERCENT SOLID

Supervisor: AMANDEEP

Analyst: jignesh

Date: 3/16/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 16:10

In Date: 03/15/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:10

Out Date: 03/16/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB101398

| <u>Lab ID</u> | <u>Client Sample ID</u> | <u>Dish#</u> | <u>Dish Wt(g)</u><br><u>(A)</u> | <u>Dish +</u><br><u>Sample Wt(g)</u><br><u>(B)</u> | <u>Dish + Dry</u><br><u>Sample Wt(g)</u><br><u>(C)</u> | <u>% Solid</u> |
|---------------|-------------------------|--------------|---------------------------------|----------------------------------------------------|--------------------------------------------------------|----------------|
| K1949-01      | BF0Y6                   | 1            | 1.13                            | 9.69                                               | 8.17                                                   | 82.2           |
| K1949-04      | BF1J5                   | 2            | 1.16                            | 9.87                                               | 7.91                                                   | 77.5           |
| K1949-06      | VHBLK02                 | 3            | 1.00                            | 2.00                                               | 2.00                                                   | 100.0          |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$