



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/20/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:10
In Date: 03/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 03/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101444

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1986-01	BF0A8	1	1.13	9.52	7.53	76.3
K1986-02	BF0A9	2	1.16	9.69	6.65	64.4
K1986-03	BF0B0	3	1.15	9.59	5.26	48.7
K1986-04	BF0B1	4	1.18	9.86	6.68	63.4
K1986-05	BF0B2	5	1.16	9.72	7.9	78.7
K1986-06	BF0B3	6	1.19	9.97	7.68	73.9
K1986-07	BF0W7	7	1.12	9.59	7.27	72.6
K1986-08	BF0Z6	8	1.11	9.71	8.34	84.1
K1986-09	BF0Z7	9	1.14	9.95	7.98	77.6
K1986-10	BF118	10	1.12	9.66	7.22	71.4
K1986-11	BF119	11	1.15	9.96	7.48	71.9
K1986-12	BF120	12	1.13	9.54	7.11	71.1
K1986-13	BF121	13	1.13	9.71	7.26	71.4
K1986-14	BF122	14	1.1	9.72	7.53	74.6
K1986-15	BF123	15	1.14	9.58	8.00	81.3
K1986-16	BF124	16	1.19	9.68	7.43	73.5
K1986-17	BF125	17	1.15	9.65	8.18	82.7
K1986-18	BF126	18	1.17	9.93	8.11	79.2
K1986-19	BF127	19	1.15	9.51	7.56	76.7
K1986-22	VHBLK02	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$