



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/21/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 03/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101480

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1988-04	BF0B4	1	1.14	9.73	7.83	77.9
K1988-05	BF0B5	2	1.11	9.6	7.99	81.0
K1988-06	BF0B6	3	1.14	9.91	7.46	72.1
K1988-07	BF0B8	4	1.14	9.88	8.19	80.7
K1988-08	BF0B9	5	1.16	9.96	6.76	63.6
K1988-09	BF0C0	6	1.17	9.92	6.37	59.4
K1988-10	BF0C1	7	1.15	9.6	7.12	70.7
K1988-11	BF0C2	8	1.13	9.72	7.27	71.5
K1988-12	BF0C3	9	1.14	9.98	8.37	81.8
K1988-13	BF0C4	10	1.14	9.54	7.55	76.3
K1988-14	BF0C5	11	1.12	9.77	8.06	80.2
K1988-15	BF0C6	12	1.15	9.64	8.03	81.0
K1988-16	BF0C8	13	1.14	9.89	8.23	81.0
K1988-17	BF0C9	14	1.11	9.64	6.21	59.8
K1988-18	BF0D0	15	1.11	9.86	5.76	53.1
K1988-19	BF0D1	16	1.15	9.72	7.52	74.3
K1988-20	BF0D2	17	1.13	9.83	8.18	81.0
K1988-21	BF0D3	18	1.13	9.89	7.6	73.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$