



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/21/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:35
In Date: 03/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101485

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2009-01	MH0YE1	1	1.00	2.00	2.00	100.0
K2009-02	MH0YE1D	2	1.00	2.00	2.00	100.0
K2009-03	MH0YE1S	3	1.00	2.00	2.00	100.0
K2009-04	MH0YS9	4	1.00	2.00	2.00	100.0
K2009-05	MH0YT0	5	1.00	2.00	2.00	100.0
K2009-06	MH0YT1	6	1.00	2.00	2.00	100.0
K2009-07	MH0YT2	7	1.00	2.00	2.00	100.0
K2009-08	MH0YT8	8	1.00	2.00	2.00	100.0
K2009-09	MH0YW0	9	1.00	2.00	2.00	100.0
K2009-10	MH0YW3	10	1.00	2.00	2.00	100.0
K2009-11	MH0YW4	11	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$