



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101511

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2011-01	20190142-COMPOSITE-A	1	1.15	9.95	5.07	44.5
K2013-01	PST-026-B123-031819	2	1.11	9.98	9.29	92.2
K2013-02	PST-026-B119-031819	3	1.11	9.57	8.57	88.2
K2013-03	PST-026-B021-031819	4	1.14	9.62	8.6	88.0
K2013-04	PST-026-B117-031819	5	1.16	9.98	8.99	88.8
K2013-05	PST-026-B120-031819	6	1.13	9.89	8.97	89.5
K2013-06	PST-026-B022-031819	7	1.15	9.64	8.96	92.0
K2013-07	PST-026-B027-031819	8	1.11	9.71	8.99	91.6
K2013-08	PST-026-B023-031819	9	1.12	9.84	9.06	91.1
K2013-09	PST-026-B024-031819	10	1.11	9.57	8.58	88.3
K2013-10	PST-026-B118-031819	11	1.15	9.7	8.95	91.2
K2013-11	PST-026-B122-031819	12	1.11	9.76	8.83	89.2
K2013-12	PST-026-SW070-031819	13	1.16	9.96	9.02	89.3
K2013-13	PST-026-SW032-031819	14	1.14	9.65	8.67	88.5
K2013-14	PST-026-SW031-031819	15	1.18	9.6	8.73	89.7
K2013-15	PST-026-SW033-031819	16	1.15	9.98	9.51	94.7
K2013-16	PST-026-SW034-031819	17	1.13	9.93	9.09	90.5
K2013-17	PST-026-SW071-031819	18	1.13	9.8	9.06	91.5
K2013-18	PST-026-SW072-031819	19	1.12	9.69	8.93	91.1
K2013-19	PST-026-SW035-031819	20	1.18	9.83	9.11	91.7
K2013-20	PST-026-SW036-031819	21	1.12	9.62	8.83	90.7
K2014-01	PST-026-SW075-031819	22	1.14	9.98	9.2	91.2
K2014-02	PST-026-SW076-031819	23	1.14	9.59	8.3	84.7
K2014-03	PST-026-SW037-031819	24	1.16	9.78	8.89	89.7
K2014-04	PST-026-SW073-031819	25	1.12	9.92	9.14	91.1
K2014-05	PST-026-SW038-031919	26	1.1	9.53	8.72	90.4
K2014-06	PST-026-SW039-031919	27	1.14	9.63	8.71	89.2
K2014-07	PST-026-SW040-031919	28	1.14	9.62	8.68	88.9
K2014-08	PST-026-SW074-031919	29	1.15	9.73	8.97	91.1
K2014-09	PST-026-SW041-031919	30	1.15	9.97	9.05	89.6
K2014-10	PST-026-SW082-031919	31	1.16	9.93	8.99	89.3
K2014-11	PST-026-B124-031919	32	1.14	9.96	9.09	90.1
K2014-12	PST-026-B121-031919	33	1.14	9.7	8.69	88.2
K2014-13	PST-026-B025-031919	34	1.15	9.68	8.7	88.5
K2014-14	PST-026-B026-031919	35	1.15	9.78	8.82	88.9
K2014-15	PST-026-SW042-031919	36	1.12	9.98	9.22	91.4
K2014-16	PST-026-SW043-031919	37	1.12	9.81	9.1	91.8
K2014-17	PST-026-SW044-031919	38	1.14	9.5	8.42	87.1
K2014-18	PST-026-SW045-031919	39	1.13	9.6	8.33	85.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101511

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2014-19	PST-026-SW077-031919	40	1.11	9.51	8.52	88.2
K2014-20	PST-026-SW078-031919	41	1.12	9.82	8.67	86.8
K2015-01	EXT-026-SW041-031919	42	1.14	9.83	8.7	87.0
K2015-02	EXT-026-SW040-031919	43	1.16	9.88	8.55	84.7
K2015-03	PST-072-SW090-031919	44	1.12	9.98	8.73	85.9
K2015-04	PST-072-B155-031919	45	1.11	9.9	9.11	91.0
K2015-06	PST-072-B157-031919	46	1.14	9.95	8.87	87.7
K2015-07	EXT-025-SW034-031919	47	1.13	9.59	8.48	86.9
K2015-08	EXT-025-SW035-031919	48	1.13	9.85	8.87	88.8
K2015-09	EXT-025-SW036-031919	49	1.14	9.64	8.87	90.9
K2015-10	EXT-025-SW037-031919	50	1.11	9.96	8.74	86.2
K2015-11	EXT-025-SW038-031919	51	1.13	9.93	8.88	88.1
K2015-12	EXT-025-SW039-031919	52	1.13	9.91	9.22	92.1
K2015-13	PST-025-B028-031919	53	1.11	9.76	9.13	92.7
K2015-14	PST-025-B029-031919	54	1.13	9.64	8.85	90.7
K2015-15	PST-025-B030-031919	55	1.13	9.8	9.14	92.4
K2015-17	PST-073-B158-031919	56	1.13	9.92	8.78	87.0
K2015-18	PST-073-B159-031919	57	1.13	9.76	8.4	84.2
K2015-19	PST-073-B160-031919	58	1.13	9.96	9.07	89.9
K2015-20	PST-074-B161-031919	59	1.12	9.58	8.83	91.1
K2040-01	S4B(2-10)	60	1.14	9.8	9.46	96.1
K2040-03	S4B(12-20)	61	1.15	9.91	8.98	89.4
K2040-05	S4B(22-28)	62	1.12	9.65	8.59	87.6
K2040-07	S4B(10)	63	1.12	9.73	8.74	88.5
K2040-09	S4B(20)	64	1.12	9.8	8.91	89.7
K2040-11	S4B(28)	65	1.14	9.77	8.49	85.2
K2043-01	WPO22SB442-190319-(12-14)	66	1.11	9.97	8.27	80.8
K2043-02	WPO22SB442-190319-FD-(12-14)	67	1.12	9.77	7.78	77.0
K2043-03	WPO22SB442-190319-(14-16)	68	1.18	9.97	8.53	83.6
K2043-04	WPO22SB429-190319-(20-22)	69	1.15	9.95	8.64	85.1
K2043-05	WPO22SB429-190319-FD-(20-22)	70	1.13	9.71	8.51	86.0
K2043-06	WPO22SB429-190319-(23-25)	71	1.16	9.77	8.19	81.6
K2043-07	WPO22SB441-190319-(23-25)	72	1.14	9.88	8.29	81.8
K2043-08	K2043-07MS	73	1.14	9.88	8.29	81.8
K2043-09	K2043-07MSD	74	1.14	9.88	8.29	81.8
K2043-10	WPO22-SB430-190319-(20-22)	75	1.17	9.94	8.77	86.7
K2043-11	WPO22SB430-190319-(23-25)	76	1.16	9.83	8.61	85.9
K2043-12	WPO22SB443-190319-(16-18)	77	1.12	9.63	8.2	83.2
K2043-13	WPO22SB443-190319-(18-20)	78	1.16	9.92	8.36	82.2



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101511

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2043-14	WPO22SB433-190319-(23-2)	79	1.16	9.5	8.00	82.0
K2043-15	WPO22SB433-190319-FD-(2)	80	1.17	9.97	8.57	84.1
K2043-16	WPO22SB433-190319-(25-2)	81	1.16	9.61	8.62	88.3
K2043-17	FIELD-QC-190319-TB2	82	1.00	2.00	2.00	100.0
K2045-01	WP022SB439-190318-(13-1)	83	1.12	9.79	8.11	80.6
K2045-02	K2045-01MS	84	1.12	9.79	8.11	80.6
K2045-03	K2045-01MSD	85	1.12	9.79	8.11	80.6
K2045-04	WP022SB439-190318-(21-2)	86	1.13	9.64	8.26	83.8
K2045-05	WP022SB439-190318-(23-2)	87	1.17	9.5	7.82	79.8
K2045-06	WP022SB439-190318-(25-2)	88	1.16	9.96	8.48	83.2
K2045-07	WP022SB438-190318-(19-2)	89	1.13	9.76	8.49	85.3
K2045-08	WP022SB438-190318-(23-2)	90	1.18	9.98	8.73	85.8
K2045-09	WP022SB438-190318-(25-2)	91	1.11	9.82	8.52	85.1
K2045-10	WP022SB435-190318-(20-2)	92	1.16	9.9	8.13	79.7
K2045-11	WP022SB434-190320-(23-2)	93	1.16	9.51	8.03	82.3
K2045-12	WP022SB434-190320-(25-2)	94	1.16	9.93	8.36	82.1
K2045-13	WP022SB436-190320-(20-2)	95	1.14	9.5	7.45	75.5
K2045-14	WP022SB436-190320-FD-(2)	96	1.17	9.97	7.87	76.1
K2045-15	WP022SB436-190320-(22-2)	97	1.14	9.52	7.76	79.0
K2045-16	WP022SB436-190320-(24-2)	98	1.17	9.97	9.04	89.4
K2045-17	WP022SB444-190320-(14-1)	99	1.1	9.5	8.41	87.0
K2045-18	WP022SB444-190320-(16-1)	100	1.1	9.91	8.4	82.9
K2045-19	WP022SB444-190320-(23-2)	101	1.13	9.56	8.2	83.9
K2045-20	FIELDQC-190320-TB3	102	1.00	2.00	2.00	100.0
K2047-01	WP022SB432-190320-(26-2)	103	1.12	9.61	8.08	82.0
K2047-02	WP022SB432-190320-(28-3)	104	1.17	9.96	8.76	86.3
K2047-03	WP022SB431-190320-(21-2)	105	1.19	9.97	8.71	85.6
K2047-04	K2047-03MS	106	1.19	9.97	8.71	85.6
K2047-05	K2047-03MSD	107	1.19	9.97	8.71	85.6
K2047-06	WP022SB431-190320-(23-2)	108	1.14	9.94	8.55	84.2
K2047-07	WP022SB445-190320-(17-1)	109	1.18	9.67	8.02	80.6
K2047-08	WP022SB445-190320-(17-1)	110	1.15	9.93	8.78	86.9
K2047-09	WP022SB445-190320-(23-2)	111	1.17	9.73	8.09	80.8
K2047-10	WP022SB434-190320-(16-1)	112	1.14	9.98	8.98	88.7
K2047-11	WP022SB435-190318-(22-2)	113	1.13	9.59	9.33	96.9
K2047-12	WP022SB435-190318-(24-2)	114	1.19	9.98	8.14	79.1
K2047-13	WP022SB435-190318-(26-2)	115	1.14	9.96	7.96	77.3
K2047-14	WP022SB437-190318-(21-2)	116	1.1	9.55	8.71	90.1
K2047-15	WP022SB437-190318-(23-2)	117	1.12	9.52	8.22	84.5



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:35
Out Date: 03/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101511

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2047-16	WP022SB440-190318-(19-2)	118	1.13	9.97	8.16	79.5
K2047-17	WP022SB440-190318-(23-2)	119	1.13	9.79	8.29	82.7
K2047-18	FIELDQC-190318-TB1	120	1.00	2.00	2.00	100.0
K2050-01	CS-5	121	1.1	9.8	8.63	86.6
K2050-02	CS-4	122	1.15	9.95	8.58	84.4
K2050-03	CS-3	123	1.14	9.92	8.24	80.9
K2050-04	CS-2	124	1.18	9.89	9.19	92.0
K2050-05	CS-1	125	1.16	9.92	8.96	89.0
K2055-01	1050	126	1.17	9.56	8.75	90.3
K2055-02	1052	127	1.13	9.63	8.7	89.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$