



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:35
In Date: 03/21/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 03/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101519

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2044-01	BF5T1	1	1.13	9.78	9.03	91.3
K2044-02	BF5T2	2	1.12	9.84	8.13	80.4
K2044-03	BF5T3	3	1.12	9.81	8.53	85.3
K2044-04	BF5T4	4	1.13	9.84	9.08	91.3
K2044-05	BF5T5	5	1.11	9.96	7.65	73.9
K2044-06	BF5T6	6	1.14	9.78	8.43	84.4
K2044-07	BF5T7	7	1.13	9.84	8.2	81.2
K2044-08	BF5T8	8	1.13	9.98	8.28	80.8
K2044-09	BF5W5	9	1.14	9.79	8.79	88.4
K2044-10	BF5W6	10	1.14	9.95	8.62	84.9
K2044-11	BF5W6MS	11	1.14	9.95	8.62	84.9
K2044-12	BF5W6MSD	12	1.14	9.95	8.62	84.9
K2044-13	BF5W2	13	1.12	9.84	9.28	93.6
K2044-14	BF5W3	14	1.11	9.82	9.11	91.8
K2044-15	BF5W4	15	1.14	9.79	6.87	66.2
K2044-16	BF5W7	16	1.14	9.55	8.36	85.9
K2044-17	VHBLK01	17	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$