



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:20
In Date: 03/22/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 03/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101577

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2077-01	MH0ZA1	1	1.00	2.00	2.00	100.0
K2077-02	MH0ZA2	2	1.00	2.00	2.00	100.0
K2077-03	MH0ZA3	3	1.00	2.00	2.00	100.0
K2077-04	MH0ZA4	4	1.00	2.00	2.00	100.0
K2077-05	MH0ZA5	5	1.00	2.00	2.00	100.0
K2077-06	MH0ZA5D	6	1.00	2.00	2.00	100.0
K2077-07	MH0ZA5S	7	1.00	2.00	2.00	100.0
K2077-08	MH0ZA6	8	1.00	2.00	2.00	100.0
K2077-09	MH0ZA7	9	1.00	2.00	2.00	100.0
K2077-10	MH0ZA8	10	1.00	2.00	2.00	100.0
K2077-11	MH0ZA9	11	1.00	2.00	2.00	100.0
K2077-12	MH0ZB0	12	1.00	2.00	2.00	100.0
K2077-13	MH0ZB1	13	1.00	2.00	2.00	100.0
K2077-14	MH0ZB2	14	1.00	2.00	2.00	100.0
K2077-15	MH0ZB3	15	1.00	2.00	2.00	100.0
K2077-16	MH0ZB4	16	1.00	2.00	2.00	100.0
K2077-17	MH0ZB5	17	1.00	2.00	2.00	100.0
K2077-18	MH0ZB6	18	1.00	2.00	2.00	100.0
K2077-19	MH0ZB7	19	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$