

TOTAL VOLATILE SOLIDS - 160.4

TEMP1 IN: 103 03/21/2019 15:20 **TEMP OUT:** 105 03/22/2019 07:00
TEMP2 IN: 105 03/22/2019 08:30 **TEMP OUT:** 105 03/22/2019 10:00
TEMP3 IN: 550 03/22/2019 11:30 **TEMP OUT:** 550 03/22/2019 13:30
TEMP4 IN: 550 03/22/2019 15:00 **TEMP OUT:** 550 03/22/2019 17:00

Run Number: LB101587

SUPERVISOR: apatel

ANALYST: Aparna

BalanceID: WC SC-3

OvenID: WC-OVEN-2

Dish #	Lab ID	Empty Dish Weight (g)	Final Empty Dish Weight (g)	Empty Dish + Sample Weight (g)	1st Dish + SampleWt Drying @103-@105°C (g)	Final Dish + SampleWt Drying @103-@105°C (g)	Dish + Samplewt Drying @550 (±50) °C (g)	Final Dish + Samplewt Drying @550 (±50) °C (g)	Weight Diff (g)	Result %
1	LB101587BL	59.6532	59.6532	59.6532	59.6532	59.6532	59.6532	59.6532	0	0
2	K2002-01	60.8424	60.8424	83.215	81.222	81.222	81.1009	81.1009	0.1211	0.6
3	K2002-01DUP	59.9835	59.9835	83.154	81.982	81.982	81.8333	81.8333	0.1487	0.7
4	K2002-03	60.8055	60.8055	83.655	82.0679	82.0679	81.112	81.112	0.9559	4.5

A = Sample Volume (ml)

B = Final Dish + Samplewt Drying @550 (±50) °C (g)

C = Final Dish + SampleWt Drying @103-@105°C (g)

D = Weight (g)

E = Final Empty Dish Weight (g)

F = Final Dish + SampleWt Drying @103-@105°C (g)

Weight D = C - B

Result mg/L = $\frac{D}{F - E} \star 100$

WORKLIST(Hardcopy Internal Chain)

WorkList Name : TVS032219

WorkList ID : 123884

Date : 03-22-2019 07:26:07

JB/0/587

Due Date	Matrix	Sample	Test	Preservative	Customer	Storage Location	Customer Sample	Collect Date	Method
03/25/2019	Solid	K2002-01	TVS	Cool 4 deg C	SKAN04	F41	S-1(2-10)-COMP	03/18/2019	160.4
03/25/2019	Solid	K2002-03	TVS	Cool 4 deg C	SKAN04	F41	S-2(5-30)-COMP	03/18/2019	160.4

Date/Time

Received by:

Relinquished by:

03/22/19, 15:00



Date/Time

Received by:

Relinquished by:

03/22/19, 15:40

