



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:30
In Date: 03/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 03/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101600

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2084-01	BF7T1	1	1.11	9.66	8.03	80.9
K2084-02	BF7T2	2	1.12	9.67	7.99	80.4
K2084-03	BF7T3	3	1.17	9.66	7.94	79.7
K2084-04	BF7W6	4	1.13	9.95	8.12	79.3
K2084-05	BF7W7	5	1.14	9.68	8.07	81.1
K2084-06	BF7W8	6	1.13	9.62	7.37	73.5
K2084-07	BF7X5	7	1.12	9.51	8.04	82.5
K2084-08	BF7X6	8	1.15	9.8	8.16	81.0
K2084-09	BF7W0	9	1.13	9.89	6.9	65.9
K2084-13	BF7X3	10	1.16	9.98	8.24	80.3
K2084-14	BF7X3MS	11	1.16	9.98	8.24	80.3
K2084-15	BF7X3MSD	12	1.16	9.98	8.24	80.3
K2084-16	BF7X4	13	1.13	9.67	8.12	81.9
K2084-19	BF7W1	14	1.11	9.63	7.94	80.2
K2084-20	BF7W2	15	1.13	9.9	8.05	78.9
K2084-21	BF7W3	16	1.14	9.83	8.33	82.7
K2084-22	BF7W4	17	1.11	9.97	8.07	78.6
K2084-23	VHBLK01	18	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$