



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:00
In Date: 03/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101637

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2036-09	WC-20190319	1	1.00	2.00	2.00	100.0
K2086-01	CS-11	2	1.00	2.00	2.00	100.0
K2086-02	CS-12	3	1.00	2.00	2.00	100.0
K2086-03	CS-13	4	1.00	2.00	2.00	100.0
K2086-04	CS-14	5	1.00	2.00	2.00	100.0
K2086-05	CS-15	6	1.00	2.00	2.00	100.0
K2086-06	CS-16	7	1.00	2.00	2.00	100.0
K2086-07	CS-18	8	1.00	2.00	2.00	100.0
K2090-01	P-5 (0-10)	9	1.12	9.59	7.97	80.9
K2090-02	P-5 (0-10)	10	1.13	9.72	7.63	75.7
K2093-01	30035-1	11	1.18	9.98	8.14	79.1
K2097-01	CS-20	12	1.18	9.74	8.48	85.3
K2097-02	CS-21	13	1.14	9.52	7.95	81.3
K2097-03	CS-22	14	1.15	9.66	8.58	87.3
K2097-04	CS-19	15	1.16	9.64	8.66	88.4
K2097-05	CS-17	16	1.18	9.83	8.61	85.9
K2097-06	CS-6	17	1.14	9.65	9.39	96.9
K2102-01	TT-001-IDWSO-20190325	18	1.18	9.97	9.35	92.9
K2102-02	TT-002-IDWSO-20190325	19	1.18	9.97	9.35	92.9
K2104-01	DRUM1-2	20	1.00	2.00	2.00	100.0
K2105-01	NJ-WATER	21	1.11	9.86	8.58	85.4
K2105-02	NJ-WATER	22	1.14	9.96	8.82	87.1
K2107-01	TP-1	23	1.15	9.95	8.61	84.8
K2107-02	TP-1	24	1.15	9.95	8.61	84.8
K2107-03	TP-2	25	1.12	9.64	8.48	86.4
K2107-04	TP-2	26	1.12	9.64	8.48	86.4
K2107-05	TP-3	27	1.17	9.5	8.82	91.8
K2107-06	TP-3	28	1.17	9.5	8.82	91.8
K2107-07	TP-4	29	1.15	9.72	9.06	92.3
K2107-08	TP-4	30	1.15	9.72	9.06	92.3
K2107-09	TP-5	31	1.16	9.54	9.02	93.8
K2107-10	TP-5	32	1.16	9.54	9.02	93.8
K2107-11	TP-6	33	1.11	9.6	8.58	88.0
K2107-12	TP-6	34	1.11	9.6	8.58	88.0
K2107-13	TP-7	35	1.16	9.94	9.00	89.3
K2107-14	TP-7	36	1.16	9.94	9.00	89.3
K2107-15	TP-8	37	1.14	9.81	9.41	95.4
K2107-16	TP-8	38	1.14	9.81	9.41	95.4
K2107-17	TP-9	39	1.18	9.97	9.28	92.2



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 3/27/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:00
In Date: 03/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 03/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101637

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>
K2107-18	TP-9	40	1.18	9.97	9.28	92.2
K2110-01	30036	41	1.1	9.54	9.43	98.7

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$