



PERCENT SOLID

Supervisor: MMDadoda

Analyst: JIGNESH

Date: 3/31/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 15:15

In Date: 03/29/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 07:28

Out Date: 03/30/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB101730

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2166-01	BF103	1	1.12	9.55	7.09	70.8
K2166-02	BF104	2	1.14	9.83	7.72	75.7
K2166-03	BF105	3	1.19	9.77	7.83	77.4
K2166-04	BF106	4	1.13	9.61	7.53	75.5
K2166-05	BF107	5	1.15	9.95	7.76	75.1
K2166-08	BF0K7	6	1.17	9.53	8.78	91.0
K2166-09	BF0L4	7	1.13	9.59	7.7	77.7
K2166-10	BF0L5	8	1.16	9.97	7.3	69.7
K2166-11	BF0L6	9	1.18	9.84	7.04	67.7
K2166-12	BF0L7	10	1.16	9.96	7.41	71.0
K2166-13	BF0M3	11	1.12	9.55	7.05	70.3
K2166-14	BF0M4	12	1.15	9.95	7.44	71.5
K2166-15	BF0M5	13	1.14	9.98	7.9	76.5
K2166-16	BF0M6	14	1.16	9.82	7.12	68.8
K2166-17	BF0N4	15	1.18	9.63	7.21	71.4
K2166-21	VHBLK01	16	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$