



PERCENT SOLID

Supervisor: MMDadoda

Analyst: JIGNESH

Date: 3/31/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 18:30

In Date: 03/29/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 03/30/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB101761

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1689-03	NB-08-032919-A	1	1.14	9.51	8.98	93.7
K1689-05	NB-08-032919-B	2	1.13	9.61	8.91	91.7
K1692-05	OK-03-032919	3	1.16	9.83	9.3	93.9
K2056-01	CS-7	4	1.00	2.00	2.00	100.0
K2056-02	CS-8	5	1.00	2.00	2.00	100.0
K2056-03	CS-9	6	1.00	2.00	2.00	100.0
K2056-04	CS-10	7	1.00	2.00	2.00	100.0
K2163-05	SVOC-GPC-BLANK	8	1.00	2.00	2.00	100.0
K2163-06	PEST-GPC-BLANK	9	1.00	2.00	2.00	100.0
K2163-07	PEST-GPC-BLANK-SPIKE	10	1.00	2.00	2.00	100.0
K2163-08	PCB-GPC-BLANK	11	1.00	2.00	2.00	100.0
K2163-09	PCB-GPC-BLANK-SPIKE	12	1.00	2.00	2.00	100.0
K2163-10	SVOC-GPC2-BLANK	13	1.00	2.00	2.00	100.0
K2163-11	PEST-GPC2-BLANK	14	1.00	2.00	2.00	100.0
K2163-12	PEST-GPC2-BLANK-SPIKE	15	1.00	2.00	2.00	100.0
K2163-13	PCB-GPC2-BLANK	16	1.00	2.00	2.00	100.0
K2163-14	PCB-GPC2-BLANK-SPIKE	17	1.00	2.00	2.00	100.0
K2171-01	GW-COMP	18	1.15	9.66	7.88	79.1
K2171-02	GW-COMP	19	1.15	9.66	7.88	79.1
K2172-01	S1	20	1.12	9.56	9.35	97.5
K2172-02	S2-SURFACE	21	1.13	9.79	9.74	99.4
K2172-03	S2-SUB-SURFACE	22	1.17	9.97	9.78	97.8
K2174-01	RBR200030	23	1.11	9.66	8.87	90.8
K2179-01	TP-4S1	24	1.15	9.68	8.84	90.2
K2179-02	TP-4S2	25	1.1	9.9	8.87	88.3
K2179-03	TP-4S	26	1.13	9.98	9.44	93.9
K2179-05	TP-4D1	27	1.12	9.67	8.42	85.4
K2179-06	TP-4D2	28	1.11	9.96	8.82	87.1
K2179-07	TP-4D	29	1.12	9.88	8.71	86.6
K2179-09	TP-8S1	30	1.12	9.68	9.14	93.7
K2179-10	TP-8S2	31	1.14	9.64	9.09	93.5
K2179-11	TP-8S	32	1.13	9.53	8.91	92.6
K2179-13	TP-8D1	33	1.13	9.5	8.88	92.6
K2179-14	TP-8D2	34	1.12	9.75	9.15	93.0
K2179-15	TP-8D	35	1.16	9.95	9.37	93.4
K2179-17	TP-5S1	36	1.11	9.92	9.02	89.8
K2179-18	TP-5S2	37	1.1	9.63	8.94	91.9
K2179-19	TP-5S	38	1.1	9.98	9.35	92.9
K2180-01	TP-7-1	39	1.14	9.58	8.86	91.5



PERCENT SOLID

Supervisor: MMDadoda

Analyst: JIGNESH

Date: 3/31/2019

OVENTEMP IN Celsius(°C): 108

Time IN: 18:30

In Date: 03/29/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103

Time OUT: 08:30

Out Date: 03/30/2019

Weight Check 1.0g: 1.00

Weight Check 10g: 10.00

BalanceID: M SC-1

Thermometer ID: %SOLIDS-OVEN

QC:LB101761

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2180-02	TP-7-2	40	1.1	9.81	9.16	92.5
K2180-03	TP-7	41	1.15	9.98	9.25	91.7
K2180-05	TP-6S1	42	1.14	9.82	9.16	92.4
K2180-06	TP-6S2	43	1.15	9.93	9.27	92.5
K2180-07	TP-6S	44	1.15	9.81	9.16	92.5
K2180-09	TP-6D1	45	1.13	9.52	8.74	90.7
K2180-10	TP-6D2	46	1.16	9.58	8.93	92.3
K2180-11	TP-6D	47	1.14	9.98	9.44	93.9
K2180-13	TP-9-1	48	1.12	9.97	9.26	92.0
K2180-14	TP-9-2	49	1.13	9.82	8.88	89.2
K2180-15	TP-9	50	1.12	9.82	9.03	90.9
K2180-17	TP-10S1	51	1.14	9.97	9.39	93.4
K2180-18	TP-10S2	52	1.13	9.84	9.32	94.0
K2180-19	TP-10	53	1.12	9.85	9.22	92.8
K2180-21	TP-10D1	54	1.13	9.62	8.81	90.5
K2180-22	TP-10D2	55	1.14	9.98	9.48	94.3
K2180-23	TP-10D	56	1.12	9.52	8.67	89.9
K2180-25	TP-1-1	57	1.12	9.99	9.38	93.1
K2180-26	TP-1-2	58	1.14	9.75	8.81	89.1
K2180-27	TP-1	59	1.13	9.84	8.86	88.7
K2180-29	TP-3-1	60	1.13	9.57	8.85	91.5
K2180-30	TP-3-2	61	1.13	9.76	8.74	88.2
K2180-31	TP-3	62	1.12	9.73	8.83	89.5

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$