



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/8/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 04/05/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 04/06/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB101900

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2251-01	MH0ZE1	1	1.00	2.00	2.00	100.0
K2251-02	MH0ZE2	2	1.00	2.00	2.00	100.0
K2251-03	MH0ZE3	3	1.00	2.00	2.00	100.0
K2251-04	MH0ZE4	4	1.00	2.00	2.00	100.0
K2251-05	MH0ZE5	5	1.00	2.00	2.00	100.0
K2251-06	MH0ZE5D	6	1.00	2.00	2.00	100.0
K2251-07	MH0ZE5S	7	1.00	2.00	2.00	100.0
K2251-08	MH0ZE6	8	1.00	2.00	2.00	100.0
K2251-09	MH0ZE7	9	1.00	2.00	2.00	100.0
K2251-10	MH0ZE8	10	1.00	2.00	2.00	100.0
K2251-11	MH0ZE9	11	1.00	2.00	2.00	100.0
K2251-12	MH0ZF0	12	1.00	2.00	2.00	100.0
K2251-13	MH0ZF1	13	1.00	2.00	2.00	100.0
K2251-14	MH0ZF2	14	1.00	2.00	2.00	100.0
K2251-15	MH0ZF3	15	1.00	2.00	2.00	100.0
K2251-16	MH0ZF4	16	1.00	2.00	2.00	100.0
K2251-17	MH0ZF5	17	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$