



PERCENT SOLID

Supervisor: apatel  
Analyst: JIGNESH  
Date: 4/8/2019

OVENTEMP IN Celsius(°C): 108  
Time IN: 16:30  
In Date: 04/05/2019  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103  
Time OUT: 08:00  
Out Date: 04/06/2019  
Weight Check 1.0g: 1.00  
Weight Check 10g: 10.00  
BalanceID: M SC-1  
Thermometer ID: %SOLIDS-OVEN

QC:LB101910

| Lab ID   | Client Sample ID | Dish# | Dish Wt(g)<br>(A) | Dish +<br>Sample Wt(g)<br>(B) | Dish + Dry<br>Sample Wt(g)<br>(C) | % Solid |
|----------|------------------|-------|-------------------|-------------------------------|-----------------------------------|---------|
| K2304-01 | BF888            | 1     | 1.16              | 9.62                          | 7.37                              | 73.4    |
| K2304-02 | BF888MS          | 2     | 1.16              | 9.62                          | 7.37                              | 73.4    |
| K2304-03 | BF888MSD         | 3     | 1.16              | 9.62                          | 7.37                              | 73.4    |
| K2304-04 | BF889            | 4     | 1.15              | 9.86                          | 7.76                              | 75.9    |
| K2304-05 | BF890            | 5     | 1.19              | 9.65                          | 8.69                              | 88.7    |
| K2304-06 | BF891            | 6     | 1.17              | 9.97                          | 7.64                              | 73.5    |
| K2304-07 | BF892            | 7     | 1.14              | 9.79                          | 7.65                              | 75.3    |
| K2304-08 | BF893            | 8     | 1.11              | 9.5                           | 8.51                              | 88.2    |
| K2304-09 | VHBLK01          | 9     | 1.00              | 2.00                          | 2.00                              | 100.0   |

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$