



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 18:10
In Date: 04/12/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:40
Out Date: 04/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102075

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2394-01	METE35	1	1.16	9.96	5.94	54.3
K2394-02	METE36	2	1.14	9.97	6.14	56.6
K2394-03	METE37	3	1.17	9.88	6.48	61.0
K2394-04	METE38	4	1.11	9.9	5.71	52.3
K2394-05	METE39	5	1.14	9.95	5.63	51.0
K2394-06	METE40	6	1.13	9.69	6.04	57.4
K2394-07	METE41	7	1.19	9.74	5.89	55.0
K2394-08	METE42	8	1.15	9.96	6.77	63.8
K2394-09	METE43	9	1.14	9.96	8.39	82.2
K2394-10	METE44	10	1.14	9.97	8.07	78.5
K2394-11	METE45	11	1.13	9.82	7.75	76.2
K2394-12	METE46	12	1.14	9.97	5.91	54.0
K2394-13	METE49	13	1.13	9.6	5.73	54.3
K2394-14	METE49D	14	1.13	9.6	5.73	54.3
K2394-15	METE49S	15	1.13	9.6	5.73	54.3
K2394-16	METE50	16	1.16	9.62	6.16	59.1
K2394-17	METE51	17	1.15	9.98	6.17	56.9
K2394-18	METE52	18	1.15	9.89	6.39	60.0
K2394-19	METE53	19	1.15	9.95	6.79	64.1
K2394-20	METE54	20	1.15	9.94	8.23	80.5
K2394-21	METE55	21	1.15	9.78	6.45	61.4
K2394-22	METE56	22	1.14	9.9	5.85	53.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$