



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 04/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102104

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2203-09	SU-01-041519-A	1	1.12	9.76	8.4	84.3
K2203-11	SU-01-041519-B	2	1.15	9.92	8.74	86.5
K2411-01	TP-16-S1	3	1.12	9.52	8.9	92.6
K2411-02	TP-16-S2	4	1.15	9.83	9.02	90.7
K2411-03	TP-16-S	5	1.17	9.66	8.94	91.5
K2411-05	TP-14-S1	6	1.15	9.9	9.41	94.4
K2411-06	TP-14-S2	7	1.15	9.96	9.58	95.7
K2411-07	TP-14-S	8	1.13	9.6	9.14	94.6
K2411-09	TP-14-D1	9	1.15	9.65	9.26	95.4
K2411-10	TP-14-D2	10	1.12	9.92	9.38	93.9
K2411-11	TP-14-D	11	1.15	9.72	9.11	92.9
K2413-01	NCC-TOPSOIL-KY	12	1.16	9.7	8.63	87.5
K2413-02	NCC-TOPSOIL-CR	13	1.11	9.56	8.5	87.5
K2414-01	FK-PS-01-044	14	1.18	9.97	8.57	84.1
K2414-02	FK-PS-01-044	15	1.16	9.63	8.68	88.8
K2417-01	FK-PS-01-045	16	1.11	9.72	7.87	78.5
K2417-02	FK-PS-01-045	17	1.19	9.74	8.54	86.0
K2422-01	TP1-S1	18	1.19	9.97	8.73	85.9
K2422-02	TP1-S2	19	1.16	9.77	8.5	85.2
K2422-03	TP1-S	20	1.14	9.94	8.68	85.7
K2422-05	TP1-D1	21	1.13	9.66	8.52	86.6
K2422-06	TP1-D2	22	1.1	9.55	8.33	85.6
K2422-07	TP1-D	23	1.14	9.91	8.74	86.7
K2422-09	TP3-S1	24	1.15	9.97	8.83	87.1
K2422-10	TP3-S2	25	1.14	9.91	8.83	87.7
K2422-11	TP3-S	26	1.15	9.85	8.31	82.3
K2423-01	TP1-A	27	1.11	9.64	7.63	76.4
K2423-02	TP1-A	28	1.11	9.64	7.63	76.4
K2423-03	TP1-B	29	1.17	9.98	7.85	75.8
K2423-04	TP1-B	30	1.17	9.98	7.85	75.8
K2423-05	TP2-A	31	1.12	9.62	7.66	76.9
K2423-06	TP2-A	32	1.12	9.62	7.66	76.9
K2423-07	TP2-B	33	1.14	9.94	8.04	78.4
K2423-08	TP2-B	34	1.14	9.94	8.04	78.4
K2423-09	TP3-A	35	1.19	9.85	8.11	79.9
K2423-10	TP3-A	36	1.19	9.85	8.11	79.9
K2423-11	TP3-B	37	1.11	9.79	8.09	80.4
K2423-12	TP3-B	38	1.11	9.79	8.09	80.4
K2423-13	TP4-A	39	1.15	9.64	8.24	83.5



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/16/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:20
In Date: 04/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:14
Out Date: 04/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102104

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2423-14	TP4-A	40	1.15	9.64	8.24	83.5
K2423-15	TP4-B	41	1.19	9.97	8.5	83.3
K2423-16	TP4-B	42	1.19	9.97	8.5	83.3
K2423-17	TP5-A	43	1.17	9.83	8.34	82.8
K2423-18	TP5-A	44	1.17	9.83	8.34	82.8
K2423-19	TP5-B	45	1.1	9.54	8.01	81.9
K2423-20	TP5-B	46	1.1	9.54	8.01	81.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$