



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/17/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 15:00
In Date: 04/16/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:32
Out Date: 04/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102117

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2402-01	GAHH3	1	1.15	9.83	8.11	80.2
K2402-02	GAHH4	2	1.16	9.88	8.55	84.7
K2402-03	GAHH5	3	1.17	9.54	8.17	83.6
K2402-04	GAHH6	4	1.16	9.7	9.18	93.9
K2402-05	GAHH7	5	1.15	9.95	8.52	83.8
K2402-06	GAHH8	6	1.17	9.91	8.83	87.6
K2402-07	GAHH9	7	1.11	9.68	9.05	92.6
K2402-08	GAHH9MS	8	1.11	9.68	9.05	92.6
K2402-09	GAHH9MSD	9	1.11	9.68	9.05	92.6
K2402-10	GAHJ0	10	1.16	9.8	8.99	90.6
K2402-11	GAHJ1	11	1.17	9.97	8.41	82.3
K2402-12	GAHJ2	12	1.18	9.62	8.58	87.7
K2402-13	GAHJ3	13	1.11	9.75	8.77	88.7
K2402-14	GAHJ4	14	1.13	9.69	9.12	93.3
K2402-15	GAHJ5	15	1.16	9.83	9.45	95.6
K2402-16	GAHJ6	16	1.17	9.74	9.47	96.8
K2402-17	GAHJ7	17	1.17	9.98	9.58	95.5
K2402-18	GAHJ8	18	1.00	2.00	2.00	100.0
K2402-19	GAHJ9	19	1.00	2.00	2.00	100.0
K2402-20	VHBLK01	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$