



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 4/19/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 04/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 04/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102202

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2403-01	A5101	1	1.15	9.96	8.93	88.3
K2403-02	A5102	2	1.14	9.71	7.3	71.9
K2403-03	A5103	3	1.19	9.78	8.82	88.8
K2403-04	A5108	4	1.00	2.00	2.00	100.0
K2403-05	A5111	5	1.00	2.00	2.00	100.0
K2403-06	A5117	6	1.00	2.00	2.00	100.0
K2403-08	A5125	7	1.11	9.77	8.55	85.9
K2403-09	A5126	8	1.15	9.6	8.69	89.2
K2403-10	A5127	9	1.14	9.85	8.64	86.1
K2403-13	VHBLK01	10	1.00	2.00	2.00	100.0
K2403-15	A5128	11	1.00	2.00	2.00	100.0
K2403-16	A5130	12	1.16	9.98	8.97	88.5
K2403-17	A5131	13	1.17	9.97	7.72	74.4
K2403-18	A5131MS	14	1.17	9.97	7.72	74.4
K2403-19	A5131MSD	15	1.17	9.97	7.72	74.4
K2403-20	A5132	16	1.14	9.62	8.61	88.1
K2403-21	A5134	17	1.12	9.65	8.8	90.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$