



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:10
In Date: 04/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 07:40
Out Date: 04/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102231

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2481-01	GAHP2	1	1.12	9.55	8.15	83.4
K2481-02	VHBLK01	2	1.00	2.00	2.00	100.0
K2481-03	GAHT8	3	1.16	9.8	7.71	75.8
K2481-04	GAHT9	4	1.16	9.97	8.2	79.9
K2481-05	GAHW0	5	1.13	9.57	7.97	81.0
K2481-06	GAHW1	6	1.11	9.71	9.08	92.7
K2481-07	GAHW2	7	1.14	9.7	8.92	90.9
K2481-08	GAHW2MS	8	1.14	9.7	8.92	90.9
K2481-09	GAHW2MSD	9	1.14	9.7	8.92	90.9
K2481-10	GAHW3	10	1.18	9.97	8.14	79.2
K2481-11	GAHW4	11	1.14	9.7	7.79	77.7
K2481-12	GAHW5	12	1.19	9.58	8.57	88.0
K2481-13	GAHW6	13	1.18	9.72	9.25	94.5
K2481-14	GAHW7	14	1.15	9.95	8.92	88.3
K2481-15	GAHW8	15	1.1	9.62	9.32	96.5
K2481-16	GAHW9	16	1.14	9.6	8.4	85.8
K2481-17	GAHX0	17	1.13	9.73	8.91	90.5
K2481-18	GAHX1	18	1.17	9.79	9.21	93.3
K2481-19	GAHX2	19	1.18	9.97	8.13	79.1
K2481-20	GAHX3	20	1.18	9.62	8.03	81.2
K2481-21	GAHX4	21	1.13	9.75	9.35	95.4
K2481-22	GAHX5	22	1.13	9.58	8.74	90.1
K2481-23	GAHX6	23	1.11	9.89	8.82	87.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$