



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 04/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:29
Out Date: 04/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102244

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2464-01	M-COMP-4.18.19	1	1.12	9.66	8.64	88.1
K2465-01	Q-COMP-4.18.19	2	1.19	9.97	9.48	94.4
K2488-05	SVOC-GPC-BLANK	3	1.00	2.00	2.00	100.0
K2488-06	PEST-GPC-BLANK	4	1.00	2.00	2.00	100.0
K2488-07	PEST-GPC-BLANK-SPIKE	5	1.00	2.00	2.00	100.0
K2488-08	PCB-GPC-BLANK	6	1.00	2.00	2.00	100.0
K2488-09	PCB-GPC-BLANK-SPIKE	7	1.00	2.00	2.00	100.0
K2488-10	SVOC-GPC2-BLANK	8	1.00	2.00	2.00	100.0
K2488-11	PEST-GPC2-BLANK	9	1.00	2.00	2.00	100.0
K2488-12	PEST-GPC2-BLANK-SPIKE	10	1.00	2.00	2.00	100.0
K2488-13	PCB-GPC2-BLANK	11	1.00	2.00	2.00	100.0
K2488-14	PCB-GPC2-BLANK -SPIKE	12	1.00	2.00	2.00	100.0
K2489-01	CONCRETE-CURB-A	13	1.00	2.00	2.00	100.0
K2489-02	CONCRETE-CURB-B	14	1.00	2.00	2.00	100.0
K2493-01	PIERSON-SOIL-01	15	1.14	9.76	6.93	67.2
K2496-01	TT-005-IDWSO-20190418	16	1.12	9.62	8.01	81.1
K2498-01	390410005	17	1.14	9.66	1.67	6.2
K2499-01	390409145	18	1.16	9.79	1.48	3.7
K2500-01	390408252	19	1.15	9.97	5.57	50.1
K2501-01	T-2	20	1.00	2.00	2.00	100.0
K2502-01	390408251	21	1.15	9.98	1.67	5.9
K2503-01	390411764	22	1.13	9.65	2.3	13.7
K2504-01	390411775	23	1.13	9.88	1.8	7.7
K2505-01	390411788	24	1.15	9.54	2.93	21.2
K2506-01	390411787	25	1.15	9.73	2.45	15.2
K2507-01	390416006	26	1.15	9.97	1.68	6.0
K2508-02	OILY-DEBRIS-220-6	27	1.00	2.00	2.00	100.0
K2508-04	OIL-220-6	28	1.00	2.00	2.00	100.0
K2508-09	OILY-DEBRIS-220-3	29	1.00	2.00	2.00	100.0
K2508-11	OIL-220-3	30	1.00	2.00	2.00	100.0
K2510-01	SOIL-PILE-A	31	1.15	9.59	8.81	90.8
K2510-02	SOIL-PILE-A	32	1.15	9.59	8.81	90.8
K2510-03	SOIL-PILE-B	33	1.13	9.55	8.42	86.6
K2510-04	SOIL-PILE-B	34	1.13	9.55	8.42	86.6
K2514-01	FK-EF-01-007-ABCD	35	1.15	9.58	8.64	88.8
K2514-02	FK-SF-01-007-ABCD	36	1.16	9.76	8.85	89.4
K2515-01	TP-15-1	37	1.16	9.55	8.55	88.1
K2515-02	TP-15-2	38	1.16	9.72	9.21	94.0
K2515-03	TP-15	39	1.14	9.78	9.17	92.9



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/22/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:40
In Date: 04/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:29
Out Date: 04/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102244

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2515-05	TP-10-1	40	1.16	9.51	9.01	94.0
K2515-06	TP-10-2	41	1.14	9.61	9.12	94.2
K2515-07	TP-10	42	1.15	9.87	9.41	94.7
K2515-09	TP-9-1	43	1.14	9.68	8.86	90.4
K2515-10	TP-9-2	44	1.15	9.9	9.08	90.6
K2515-11	TP-9	45	1.14	9.95	9.21	91.6
K2515-13	TP-12-1	46	1.15	9.91	9.06	90.3
K2515-14	TP-12-2	47	1.15	9.65	9.17	94.4
K2515-15	TP-12	48	1.14	9.95	9.48	94.7
K2515-17	TP-13-1	49	1.15	9.67	8.57	87.1
K2515-18	TP-13-2	50	1.14	9.56	8.84	91.4
K2515-19	TP-13	51	1.14	9.95	9.07	90.0
K2515-21	TP-16-1	52	1.15	9.61	8.81	90.5
K2515-22	TP-16-2	53	1.16	9.77	9.22	93.6
K2515-23	TP-16	54	1.16	9.98	9.32	92.5

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$