



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/24/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 04/23/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:10
Out Date: 04/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102304

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2200-07	EO-02-042319-A	1	1.13	9.62	8.61	88.1
K2200-09	EO-02-042319-B	2	1.16	9.73	8.65	87.4
K2200-11	EO-02-042319-C	3	1.15	9.58	8.46	86.7
K2538-03	NB-20523	4	1.00	2.00	2.00	100.0
K2539-01	FK-SF-03-006-A	5	1.14	9.68	8.74	89.0
K2539-02	FK-SF-03-006-A	6	1.15	9.51	8.61	89.2
K2539-03	FK-SF-03-006-B	7	1.16	9.82	8.81	88.3
K2539-04	FK-SF-03-006-B	8	1.17	9.73	8.84	89.6
K2540-01	FK-SF-03-007-A	9	1.14	9.6	8.65	88.8
K2540-02	FK-SF-03-007-A	10	1.16	9.75	8.86	89.6
K2540-03	FK-SF-03-007-B	11	1.12	9.56	8.64	89.1
K2540-04	FK-SF-03-007-B	12	1.14	9.62	8.76	89.9
K2542-01	TP-6-1	13	1.16	9.63	8.4	85.5
K2542-02	TP-6-2	14	1.15	9.71	8.66	87.7
K2542-03	TP-6	15	1.14	9.95	9.00	89.2
K2542-05	TP-1-1	16	1.15	9.82	8.78	88.0
K2542-06	TP-1-2	17	1.14	9.98	9.51	94.7
K2542-07	TP-1	18	1.14	9.97	9.24	91.7
K2542-09	TP-7-1	19	1.14	9.74	8.3	83.3
K2542-10	TP-7-2	20	1.14	9.71	8.79	89.3
K2542-11	TP-7	21	1.14	9.86	8.74	87.2
K2542-13	TP-11-1	22	1.14	9.8	8.76	88.0
K2542-14	TP-11-2	23	1.14	9.82	9.2	92.9
K2542-15	TP-11	24	1.13	9.67	9.08	93.1
K2542-17	TP-8-1	25	1.15	9.65	8.6	87.6
K2542-18	TP-8-2	26	1.15	9.72	8.47	85.4
K2542-19	TP-8	27	1.16	9.65	8.62	87.9

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$