



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 04/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 04/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102325

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2203-17	SU-03-042419-A	33	1.14	9.69	8.69	88.3
K2203-19	SU-03-042419-B	34	1.15	9.95	8.89	88.0
K2533-01	PST-028-B001-042219	1	1.14	9.87	9.18	92.1
K2533-02	PST-028-B010-042219	2	1.12	9.71	9.08	92.7
K2533-03	PST-028-B005-042219	3	1.18	9.58	8.73	89.9
K2533-04	PST-028-B006-042219	4	1.15	9.74	8.78	88.8
K2533-05	PST-028-B008-042219	5	1.15	9.95	8.5	83.5
K2533-06	PST-027-B109-042219	6	1.14	9.88	8.79	87.5
K2533-07	PST-072-B156-042219	7	1.11	9.56	8.66	89.3
K2533-08	PST-025-B031-042219	8	1.13	9.51	8.59	89.0
K2533-09	PST-023-B183-042219	9	1.16	9.8	8.89	89.5
K2533-10	PST-023-B184-042219	10	1.16	9.97	9.41	93.6
K2533-11	PST-023-B182-042219	11	1.18	9.58	8.89	91.8
K2533-12	PST-023-B181-042219	12	1.15	9.55	8.57	88.3
K2533-13	PST-028-SW005-042219	13	1.17	9.79	8.91	89.8
K2533-14	EXT-028-SW001-042219	14	1.1	9.85	8.76	87.5
K2533-15	EXT-028-SW023-042219	15	1.13	9.71	8.67	87.9
K2533-16	EXT-028-SW022-042219	16	1.16	9.98	8.61	84.5
K2533-17	EXT-028-SW021-042219	17	1.14	9.94	9.06	90.0
K2533-18	EXT-027-SW020-042219	18	1.11	9.96	8.89	87.9
K2533-19	EXT-027-SW019-042219	19	1.12	9.63	8.69	89.0
K2533-20	EXT-027-SW018-042219	20	1.17	9.65	8.24	83.4
K2534-01	EXT-074-SW065-042219	21	1.16	9.97	8.36	81.7
K2534-02	EXT-074-SW066-042219	22	1.15	9.63	8.14	82.4
K2534-03	EXT-074-SW063-042219	23	1.15	9.66	8.45	85.8
K2534-04	EXT-073-SW067-042219	24	1.15	9.66	8.51	86.5
K2534-05	EXT-073-SW061-042219	25	1.15	9.98	8.65	84.9
K2534-06	EXT-073-SW060-042219	26	1.14	9.56	8.46	86.9
K2534-07	EXT-073-SW059-042219	27	1.17	9.81	8.45	84.3
K2534-08	EXT-073-SW058-042219	28	1.14	9.58	8.43	86.4
K2534-09	EXT-025-SW031-042219	29	1.16	9.64	8.69	88.8
K2534-10	EXT-025-SW030-042219	30	1.15	9.87	8.89	88.8
K2544-01	135-20	35	1.13	9.96	8.72	86.0
K2552-01	S23432-01-71391	31	1.00	2.00	2.00	100.0
K2552-02	S23432-02-71392	32	1.00	2.00	2.00	100.0
K2553-01	RBR200034	36	1.1	9.59	8.31	84.9
K2553-02	RBR200034	37	1.1	9.59	8.31	84.9
K2553-04	RBR200051	38	1.18	9.9	8.36	82.3
K2553-05	RBR200051	39	1.18	9.9	8.36	82.3



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/25/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 04/24/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 04/25/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102325

<u>Lab ID</u>	<u>Client Sample ID</u>	<u>Dish#</u>	<u>Dish Wt(g)</u> <u>(A)</u>	<u>Dish +</u> <u>Sample Wt(g)</u> <u>(B)</u>	<u>Dish + Dry</u> <u>Sample Wt(g)</u> <u>(C)</u>	<u>% Solid</u>
K2554-01	RT3860	40	1.16	9.66	8.77	89.5
K2559-01	390418012	41	1.18	9.64	2.73	18.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$