



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 04/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 04/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102382

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2551-01	PT-AN-SOIL	1	1.00	2.00	2.00	100.0
K2551-02	PT-CORR-SOIL	2	1.00	2.00	2.00	100.0
K2551-03	PT-CORR-SOIL	3	1.00	2.00	2.00	100.0
K2551-04	PT-CN-SOIL	4	1.00	2.00	2.00	100.0
K2551-05	PT-CN-SOIL	5	1.00	2.00	2.00	100.0
K2551-06	PT-FP-SOIL	6	1.00	2.00	2.00	100.0
K2551-07	PT-CR6-SOIL	7	1.00	2.00	2.00	100.0
K2551-08	PT-NUT-SOIL	8	1.00	2.00	2.00	100.0
K2551-09	PT-NUT-SOIL	9	1.00	2.00	2.00	100.0
K2551-10	PT-OGR-SOIL	10	1.00	2.00	2.00	100.0
K2551-11	PT-MET-SOIL	11	1.00	2.00	2.00	100.0
K2551-12	PT-BNA-SOIL	12	1.00	2.00	2.00	100.0
K2551-13	PT-PAH-SOIL	13	1.00	2.00	2.00	100.0
K2551-14	PT-DIES-SOIL	14	1.00	2.00	2.00	100.0
K2551-15	PT-GAS-SOIL	15	1.00	2.00	2.00	100.0
K2551-16	PT-NJEPH-SOIL	16	1.00	2.00	2.00	100.0
K2551-17	PT-HERB-SOIL	17	1.00	2.00	2.00	100.0
K2551-18	PT-PCB-SOIL	18	1.00	2.00	2.00	100.0
K2551-19	PT-PEST-SOIL	19	1.00	2.00	2.00	100.0
K2551-20	PT-CHLR-SOIL	20	1.00	2.00	2.00	100.0
K2551-21	PT-TXP-SOIL	21	1.00	2.00	2.00	100.0
K2551-22	PT-VOA-SOIL	22	1.00	2.00	2.00	100.0
K2551-23	PT-SOL-SOIL	23	1.00	2.00	2.00	100.0
K2582-01	SSSI-0419-S-DH-6-T	24	1.14	9.89	9.1	91.0
K2582-02	SSSI-0419-S-DH-6-T	25	1.18	9.9	9.38	94.0
K2582-03	SSSI-0419-S-DH-6-B	26	1.15	9.7	9.01	91.9
K2582-04	SSSI-0419-S-DH-6-B	27	1.1	9.82	9.27	93.7
K2583-02	30055-23933	28	1.13	9.63	9.32	96.4
K2584-01	SSSI-0419-S-DH-7	29	1.19	9.89	8.93	89.0
K2584-02	SSSI-0419-S-DH-7	30	1.18	9.61	9.04	93.2
K2584-03	SSSI-0419-S-DH-8	31	1.15	9.51	8.57	88.8
K2584-04	SSSI-0419-S-DH-8	32	1.14	9.8	8.96	90.3
K2587-01	EXT-004-SW069-042319	33	1.13	9.85	8.62	85.9
K2587-02	PST-003-B173-042319	34	1.16	9.71	8.8	89.4
K2587-03	PST-003-B172-042319	35	1.16	9.62	8.89	91.4
K2587-04	PST-003-B171-042319	36	1.14	9.84	8.74	87.4
K2587-05	PST-003-B176-042319	37	1.15	9.88	9.06	90.6
K2587-06	PST-003-B177-042319	38	1.15	9.96	9.51	94.9
K2587-07	PST-003-B178-042319	39	1.16	9.73	9.22	94.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 04/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 04/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102382

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2587-08	PST-004-B168-042319	40	1.16	9.5	8.83	92.0
K2587-09	PST-004-B166-042319	41	1.16	9.64	8.93	91.6
K2587-10	PST-004-B165-042319	42	1.16	9.56	8.81	91.1
K2587-11	PST-004-SW092-042319	43	1.17	9.9	8.86	88.1
K2587-12	PST-004-SW022-042319	44	1.18	9.9	8.91	88.6
K2587-13	PST-004-SW023-042319	45	1.1	9.54	8.58	88.6
K2587-14	PST-004-SW091-042319	46	1.18	9.86	8.57	85.1
K2587-15	PST-004-B169-042319	47	1.18	9.53	8.53	88.0
K2587-16	PST-004-B162-042319	48	1.16	9.59	8.91	91.9
K2587-17	PST-004-B163-042319	49	1.15	9.97	9.34	92.9
K2587-18	PST-004-B164-042319	50	1.14	9.65	8.93	91.5
K2587-19	PST-004-B202-042319	51	1.17	9.7	8.88	90.4
K2587-20	EXT-004-SW068-042319	52	1.17	9.98	8.89	87.6
K2588-01	PST-004-B167-042319	53	1.16	9.81	9.07	91.4
K2588-02	PST-004-B170-042319	54	1.14	9.79	8.8	88.6
K2588-03	PST-001-B016-042519	55	1.16	9.56	8.89	92.0
K2588-04	PST-001-B193-042519	56	1.18	9.66	8.91	91.2
K2588-05	PST-001-B015-042519	57	1.14	9.82	8.9	89.4
K2588-06	PST-001-B012-042519	58	1.19	9.55	8.6	88.6
K2588-07	PST-001-B013-042519	59	1.15	9.64	9.08	93.4
K2588-08	PST-001-B192-042519	60	1.15	9.8	9.11	92.0
K2588-09	PST-001-B014-042519	61	1.17	9.97	9.47	94.3
K2588-10	PST-001-B194-042519	62	1.18	9.9	9.27	92.8
K2589-01	PST-001-B203-042519	63	1.14	9.52	8.75	90.8
K2589-02	PST-001-B195-042519	64	1.13	9.58	8.6	88.4
K2589-03	PST-001-B196-042519	65	1.14	9.69	8.61	87.4
K2589-04	PST-001-B197-042519	66	1.14	9.97	8.09	78.7
K2589-05	PST-001-SW018-042419	67	1.19	9.97	9.31	92.5
K2589-06	PST-001-SW017-042419	68	1.11	9.55	8.72	90.2
K2589-07	PST-001-SW098-042419	69	1.17	9.65	8.83	90.3
K2589-08	PST-001-SW016-042419	70	1.17	9.58	8.76	90.2
K2589-09	PST-001-B019-042419	71	1.19	9.97	9.42	93.7
K2589-10	PST-001-B201-042419	72	1.15	9.88	9.24	92.7
K2589-11	PST-001-B190-042419	73	1.18	9.64	8.99	92.3
K2589-12	PST-001-B191-042419	74	1.12	9.56	8.87	91.8
K2589-13	PST-001-B189-042419	75	1.19	9.73	8.9	90.3
K2589-14	PST-001-B188-042419	76	1.18	9.68	8.99	91.9
K2589-15	PST-001-B186-042419	77	1.13	9.57	8.88	91.8
K2589-16	PST-001-B187-042419	78	1.17	9.73	9.15	93.2



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/29/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 04/26/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:25
Out Date: 04/27/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102382

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2589-17	PST-001-B185-042419	79	1.1	9.95	9.24	92.0
K2591-05	SVOC-GPC-BLANK	80	1.00	2.00	2.00	100.0
K2591-06	PEST-GPC-BLANK	81	1.00	2.00	2.00	100.0
K2591-07	PEST-GPC-BLANK-SPIKE	82	1.00	2.00	2.00	100.0
K2591-08	PCB-GPC-BLANK	83	1.00	2.00	2.00	100.0
K2591-09	PCB-GPC-BLANK-SPIKE	84	1.00	2.00	2.00	100.0
K2591-10	SVOC-GPC2-BLANK	85	1.00	2.00	2.00	100.0
K2591-11	PEST-GPC2-BLANK	86	1.00	2.00	2.00	100.0
K2591-12	PEST-GPC2-BLANK-SPIKE	87	1.00	2.00	2.00	100.0
K2591-13	PCB-GPC2-BLANK	88	1.00	2.00	2.00	100.0
K2591-14	PCB-GPC2-BLANK -SPIKE	89	1.00	2.00	2.00	100.0
K2600-01	T-1	90	1.16	9.76	7.53	74.1
K2600-02	T-2	91	1.15	9.91	7.79	75.8
K2600-03	T-3	92	1.19	9.71	7.71	76.5
K2601-01	FK-EF-01-008-ABCD	93	1.19	9.69	8.53	86.4
K2601-02	FK-EF-01-008-ABCD	94	1.1	9.83	8.65	86.5

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$