



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 4/30/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:10
In Date: 04/29/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:11
Out Date: 04/30/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102417

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2194-01	JC-03-042919-A	1	1.15	9.85	9.05	90.8
K2194-03	JC-03-042919-B	2	1.19	9.88	9.07	90.7
K2194-05	JC-03-042919-C	3	1.18	9.81	9.08	91.5
K2194-07	JC-03-042919-D	4	1.15	9.64	8.94	91.8
K2194-09	JC-03-042919-E	5	1.15	9.98	8.98	88.7
K2194-11	JC-03-042919-F	6	1.12	9.79	9.12	92.3
K2194-13	JC-03-042919-G	7	1.19	9.72	8.68	87.8
K2194-15	JC-03-042919-H	8	1.14	9.63	8.75	89.6
K2194-17	JC-03-042919-I	9	1.15	9.95	9.26	92.2
K2199-17	OR-02-042919-A	10	1.11	9.64	8.73	89.3
K2199-19	OR-02-042919-B	11	1.13	9.82	9.23	93.2
K2608-01	DRUM-A-E	12	1.00	2.00	2.00	100.0
K2612-01	B-1	13	1.18	9.52	9.00	93.8
K2612-02	B-2	14	1.16	9.62	9.11	94.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$