



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 05/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102457

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2602-01	BFH48	1	1.13	9.66	8.54	86.9
K2602-02	BFH48MS	2	1.13	9.66	8.54	86.9
K2602-03	BFH48MSD	3	1.13	9.66	8.54	86.9
K2602-04	BFH59	4	1.16	9.7	8.7	88.3
K2602-05	BFH61	5	1.18	9.64	8.48	86.3
K2602-06	BFH69	6	1.14	9.89	8.74	86.9
K2602-07	VHBLK01	7	1.00	2.00	2.00	100.0
K2602-08	BFH52	8	1.14	9.71	8.96	91.2
K2602-09	BFH58	9	1.18	9.56	8.52	87.6
K2602-10	BFH65	10	1.17	9.73	9.33	95.3
K2602-11	BFH70	11	1.16	9.84	9.47	95.7
K2602-12	BFHD7	12	1.17	9.67	9.21	94.6
K2602-13	BFHD8	13	1.12	9.69	9.22	94.5
K2602-14	BFHF1	14	1.1	9.83	9.33	94.3
K2602-15	BFHF2	15	1.18	9.59	8.42	86.1
K2602-16	BFHK9	16	1.18	9.57	8.54	87.7

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$