



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:20
In Date: 05/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:32
Out Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102458

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2604-01	GAJ75	1	1.16	9.72	9.09	92.6
K2604-02	GAJ76	2	1.18	9.98	9.34	92.7
K2604-03	GAJ77	3	1.16	9.72	8.99	91.5
K2604-04	GAJ78	4	1.14	9.96	9.05	89.7
K2604-05	GAJ79	5	1.15	9.59	9.34	97.0
K2604-06	GAJ80	6	1.14	9.87	8.91	89.0
K2604-07	GAJ81	7	1.15	9.51	9.12	95.3
K2604-08	GAJ82	8	1.14	9.85	9.37	94.5
K2604-09	GAJ83	9	1.13	9.62	8.74	89.6
K2604-10	GAJ84	10	1.13	9.69	9.32	95.7
K2604-11	GAJ85	11	1.14	9.61	7.65	76.9
K2604-12	GAJ86	12	1.14	9.6	8.07	81.9
K2604-13	GAJ87	13	1.15	9.72	8.68	87.9
K2604-14	GAJ88	14	1.13	9.72	8.68	87.9
K2604-15	GAJ88MS	15	1.13	9.72	8.68	87.9
K2604-16	GAJ88MSD	16	1.13	9.72	8.68	87.9
K2604-17	VHBLK01	17	1.00	2.00	2.00	100.0
K2604-18	GAJ89	18	1.13	9.97	9.2	91.3
K2604-19	GAJ90	19	1.16	9.95	9.57	95.7
K2604-20	GAJ91	20	1.16	9.97	8.71	85.7
K2604-21	GAJ92	21	1.15	9.8	9.11	92.0
K2604-22	GAJ93	22	1.14	9.98	8.29	80.9
K2604-23	GAJ94	23	1.13	9.6	8.37	85.5

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$