



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:35
In Date: 05/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:37
Out Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102464

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2605-01	MA5110	1	1.00	2.00	2.00	100.0
K2605-02	MA5117	2	1.00	2.00	2.00	100.0
K2605-03	MA5138	3	1.18	9.96	8.84	87.2
K2605-04	MA5139	4	1.19	9.96	8.92	88.1
K2605-05	MA5140	5	1.14	9.98	8.93	88.1
K2605-06	MA5141	6	1.18	9.97	8.69	85.4
K2605-08	MA5144	7	1.15	9.63	9.08	93.5
K2605-09	MA5145	8	1.14	9.84	8.82	88.3
K2605-10	MA5146	9	1.13	9.96	8.74	86.2
K2605-11	MA5146D	10	1.13	9.96	8.74	86.2
K2605-12	MA5146S	11	1.13	9.96	8.74	86.2
K2605-13	MA5147	12	1.15	9.53	8.78	91.1
K2605-14	MA5148	13	1.13	9.85	8.45	83.9
K2605-15	MA5149	14	1.14	9.78	8.85	89.2
K2605-16	MA5150	15	1.15	9.66	8.67	88.4
K2605-17	MA5151	16	1.13	9.72	8.82	89.5
K2605-18	MA5152	17	1.14	9.94	8.57	84.4
K2605-19	MA5105	18	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$