



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 05/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:42
Out Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102467

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2633-01	GAJ95	1	1.15	9.51	8.51	88.0
K2633-02	GAJ95MS	2	1.15	9.51	8.51	88.0
K2633-03	GAJ95MSD	3	1.15	9.51	8.51	88.0
K2633-04	GAJ96	4	1.16	9.63	9.11	93.9
K2633-05	GAJ97	5	1.15	9.63	9.2	94.9
K2633-06	GAJ98	6	1.15	9.82	9.56	97.0
K2633-07	GAJ99	7	1.14	9.92	9.47	94.9
K2633-08	GAJA0	8	1.16	9.61	9.34	96.8
K2633-09	GAJA1	9	1.17	9.69	8.12	81.6
K2633-10	GAJA2	10	1.15	9.98	9.1	90.0
K2633-11	GAJA3	11	1.15	9.75	8.15	81.4
K2633-12	GAJA4	12	1.15	9.97	8.55	83.9
K2633-13	GAJA5	13	1.15	9.72	8.5	85.8
K2633-14	GAJA6	14	1.14	9.88	9.47	95.3
K2633-15	GAJA7	15	1.15	9.59	8.19	83.4
K2633-16	GAJA8	16	1.15	9.93	8.94	88.7
K2633-17	GAJA9	17	1.13	9.74	9.52	97.4
K2633-18	GAJB0	18	1.14	9.72	9.36	95.8
K2633-19	GAJB1	19	1.15	9.7	7.76	77.3
K2633-20	VHBLK01	20	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$