



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/2/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:50
In Date: 05/01/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:05
Out Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102486

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2627-01	WC-1	1	1.18	9.93	8.88	88.0
K2650-04	OILY-DEBRIS	2	1.00	2.00	2.00	100.0
K2652-01	DRUM-1-4	3	1.15	9.95	7.73	74.8
K2654-01	OILY-STONES	4	1.00	2.00	2.00	100.0
K2654-02	OILY-STONES-A	5	1.00	2.00	2.00	100.0
K2654-03	OILY-STONES-B	6	1.00	2.00	2.00	100.0
K2654-04	OILY-STONES	7	1.00	2.00	2.00	100.0
K2656-01	390425737	8	1.15	9.97	3.34	24.8
K2657-01	390430013	9	1.14	9.96	1.39	2.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$