



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102518

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K1411-04	ETGI-RO-294	1	1.16	9.63	8.88	91.1
K2635-01	AU-05-050219-A	2	1.13	9.84	9.05	90.9
K2635-03	AU-05-050219-B	3	1.14	9.64	8.88	91.1
K2635-05	AU-05-050219-C	4	1.14	9.8	9.00	90.8
K2635-07	AU-05-050219-D	5	1.14	9.56	8.77	90.6
K2642-01	OK-02-050219	6	1.14	9.92	9.44	94.5
K2647-01	SU-04-050219-A	7	1.13	9.97	8.83	87.1
K2647-03	SU-04-050219-B	8	1.15	9.68	8.71	88.6
K2667-01	SP-1-A	9	1.12	9.65	7.82	78.5
K2667-03	SP-1-B	10	1.12	9.65	8.48	86.3
K2668-01	TP-1	11	1.15	9.9	9.41	94.4
K2668-02	TP-1	12	1.15	9.9	9.41	94.4
K2668-03	TP-2	13	1.15	9.53	8.54	88.2
K2668-04	TP-2	14	1.15	9.53	8.54	88.2
K2668-05	TP-3	15	1.11	9.71	8.14	81.7
K2668-06	TP-3	16	1.11	9.71	8.14	81.7
K2669-01	OILY-DEBRIS-AB	17	1.00	2.00	2.00	100.0
K2669-02	OIL-DRUM	18	1.00	2.00	2.00	100.0
K2672-01	SB-01-1.0-1.5	19	1.15	9.97	9.2	91.3
K2672-02	SB-01-COMP	20	1.15	9.89	8.96	89.4
K2672-03	SB-02-1.0-1.5	21	1.16	9.7	8.87	90.3
K2672-04	SB-02-COMP	22	1.14	9.84	9.15	92.1
K2672-05	SB-03-1.0-1.5	23	1.14	9.66	7.83	78.5
K2672-06	SB-03-COMP	24	1.12	9.98	8.92	88.0
K2672-07	SB-04-1.0-1.5	25	1.14	9.75	8.4	84.3
K2672-08	SB-04-COMP	26	1.15	9.73	8.33	83.7
K2672-09	SB-05-1.5-2.0	27	1.15	9.53	7.73	78.5
K2672-10	SB-05-COMP	28	1.16	9.68	8.36	84.5
K2672-11	SB-06-1.5-2.0	29	1.14	9.79	8.1	80.5
K2672-12	SB-06-COMP	30	1.16	9.6	8.41	85.9
K2672-13	SB-07-1.5-2.0	31	1.15	9.96	8.68	85.5
K2672-14	SB-07-COMP	32	1.14	9.9	8.98	89.5
K2672-15	SB-08-COMP	33	1.15	9.97	8.59	84.4
K2672-16	SB-08-COMP	34	1.15	9.97	8.59	84.4
K2672-17	SB-09-1.0-1.5	35	1.14	9.74	8.9	90.2
K2672-18	SB-09-COMP	36	1.15	9.57	8.6	88.5
K2672-19	SB-10-1.0-1.5	37	1.14	9.53	8.63	89.3
K2672-20	SB-10-COMP	38	1.15	9.6	8.81	90.7
K2673-01	SB-11-1.5-2.0	39	1.15	9.95	9.09	90.2



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/15/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:35
In Date: 05/02/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102518

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2673-02	SB-11-COMP	40	1.14	9.56	8.8	91.0
K2673-03	SB-12-2.0-2.5	41	1.15	9.62	8.21	83.4
K2673-04	SB-12-COMP	42	1.14	9.59	8.26	84.3
K2673-05	SB-13-1.5-2.0	43	1.13	9.92	7.57	73.3
K2673-06	SB-13-COMP	44	1.15	9.61	8.21	83.5
K2673-07	SB-14-1.0-1.5	45	1.14	9.98	8.56	83.9
K2673-08	SB-14-COMP	46	1.14	9.64	8.26	83.8
K2673-09	SB-15-1.0-1.5	47	1.15	9.96	8.84	87.3
K2673-10	SB-15-COMP	48	1.15	9.96	9.09	90.1
K2673-11	SB-16-2.0-2.5	49	1.15	9.58	7.61	76.6
K2673-12	SB-16-COMP	50	1.15	9.94	8.92	88.4
K2673-13	SB-17-1.0-1.5	51	1.14	9.72	8.6	86.9
K2673-14	SB-17-COMP	52	1.15	9.65	8.55	87.1
K2673-15	SP-01-COMP	53	1.15	9.77	8.8	88.7
K2673-16	SP-02-COMP	54	1.14	9.98	9.14	90.5
K2674-01	FILTERBAG-RO	55	1.00	2.00	2.00	100.0
K2677-01	FK-EF-01-007-ABCD	56	1.16	9.72	8.62	87.1
K2677-02	FK-EF-01-007-ABCD	57	1.15	9.7	8.66	87.8
K2678-01	FK-0D015-03-003	58	1.16	9.7	9.37	96.1
K2678-02	FK-0D015-03-003	59	1.15	9.64	9.31	96.1
K2679-01	T-1	60	1.14	9.8	8.11	80.5
K2679-02	T-2	61	1.14	9.98	8.91	87.9
K2679-03	T-3	62	1.13	9.63	7.99	80.7
K2679-04	T-4	63	1.15	9.74	8.00	79.7
K2679-05	T-5	64	1.15	9.7	8.64	87.6

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$