



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/6/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 05/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 05/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102540

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2690-01	A5106	1	1.00	2.00	2.00	100.0
K2690-02	A5112	2	1.00	2.00	2.00	100.0
K2690-03	A5116	3	1.00	2.00	2.00	100.0
K2690-04	A5150	4	1.14	9.55	8.52	87.8
K2690-05	A5151	5	1.13	9.78	8.87	89.5
K2690-06	A5152	6	1.14	9.8	8.51	85.1
K2690-07	A5153	7	1.14	9.81	9.27	93.8
K2690-08	A5154	8	1.16	9.97	9.21	91.4
K2690-09	A5155	9	1.13	9.66	7.83	78.5
K2690-10	A5156	10	1.14	9.98	8.1	78.7
K2690-11	A5157	11	1.13	9.77	9.36	95.3
K2690-12	A5158	12	1.13	9.96	8.99	89.0
K2690-13	A5159	13	1.15	9.6	8.4	85.8
K2690-14	A5160	14	1.16	9.82	9.25	93.4
K2690-15	A5161	15	1.14	9.71	8.46	85.4
K2690-16	A5162	16	1.15	9.7	8.59	87.0
K2690-17	A5163	17	1.14	9.75	9.21	93.7
K2690-18	A5164	18	1.16	9.83	8.8	88.1
K2690-19	A5165	19	1.16	9.96	9.08	90.0
K2690-20	A5166	20	1.15	9.61	9.14	94.4
K2690-21	A5167	21	1.16	9.77	8.75	88.2
K2690-22	A5167MS	22	1.16	9.77	8.75	88.2
K2690-23	A5167MSD	23	1.16	9.77	8.75	88.2
K2690-24	A5168	24	1.16	9.97	8.93	88.2
K2690-27	VHBLK02	25	1.00	2.00	2.00	100.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$