



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/6/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:15
In Date: 05/03/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 05/04/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102541

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2688-01	MA5104	1	1.00	2.00	2.00	100.0
K2688-02	MA5113	2	1.00	2.00	2.00	100.0
K2688-03	MA5116	3	1.00	2.00	2.00	100.0
K2688-04	MA5153	4	1.14	9.81	9.27	93.8
K2688-05	MA5154	5	1.16	9.97	9.21	91.4
K2688-06	MA5155	6	1.13	9.66	7.83	78.5
K2688-07	MA5156	7	1.14	9.98	8.1	78.7
K2688-08	MA5157	8	1.13	9.77	9.36	95.3
K2688-09	MA5158	9	1.13	9.96	8.99	89.0
K2688-10	MA5159	10	1.15	9.6	8.4	85.8
K2688-11	MA5160	11	1.16	9.82	9.25	93.4
K2688-12	MA5161	12	1.14	9.71	8.46	85.4
K2688-13	MA5162	13	1.15	9.7	8.59	87.0
K2688-14	MA5163	14	1.14	9.75	9.21	93.7
K2688-15	MA5164	15	1.16	9.83	8.8	88.1
K2688-16	MA5165	16	1.16	9.96	9.08	90.0
K2688-17	MA5166	17	1.15	9.61	9.14	94.4
K2688-18	MA5167	18	1.16	9.77	8.75	88.2
K2688-19	MA5167D	19	1.16	9.77	8.75	88.2
K2688-20	MA5167S	20	1.16	9.77	8.75	88.2
K2688-21	MA5168	21	1.16	9.97	8.93	88.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$