



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/8/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:25
In Date: 05/07/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 05/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102594

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2717-01	P026-SS001-0002-01	1	1.15	9.96	8.65	85.1
K2717-02	P026-SS001-0206-01	2	1.16	9.98	9.07	89.7
K2717-03	P026-SS001-0612-01	3	1.17	9.68	9.02	92.2
K2717-04	P026-SS001-1218-01	4	1.16	9.53	8.64	89.4
K2717-05	P026-SS001-1824-01	5	1.16	9.76	8.75	88.3
K2717-06	P026-SS002-0002-01	6	1.16	9.65	7.25	71.7
K2717-07	P026-SS002-0206-01	7	1.16	9.61	7.7	77.4
K2717-08	P026-SS002-0612-01	8	1.14	9.78	8.04	79.9
K2717-09	P026-SS002-1218-01	9	1.15	9.59	7.92	80.2
K2717-10	K2717-09MS	10	1.15	9.59	7.92	80.2
K2717-11	K2717-09MSD	11	1.15	9.59	7.92	80.2
K2717-12	P026-SS002-1218-02	12	1.14	9.7	8.02	80.4
K2717-13	P026-SS002-1824-01	13	1.15	9.93	8.54	84.2
K2717-14	P026-SS003-0002-01	14	1.15	9.73	7.14	69.8
K2717-15	P026-SS003-0206-01	15	1.14	9.82	7.83	77.1
K2717-16	P026-SS003-0612-01	16	1.17	9.51	7.51	76.0
K2717-17	P026-SS003-1218-01	17	1.15	9.56	7.56	76.2
K2717-18	P026-SS003-1824-01	18	1.14	9.92	7.92	77.2
K2717-19	P026-SS004-0002-01	19	1.15	9.68	8.22	82.9
K2717-20	P026-SS004-0206-01	20	1.14	9.67	8.82	90.0
K2717-21	P026-SS004-0612-01	21	1.14	9.69	9.05	92.5

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$