



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:35
In Date: 05/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102604

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2720-01	MH0AQ3	1	1.14	9.63	8.46	86.2
K2720-02	MH0AQ4	2	1.15	9.8	8.23	81.8
K2720-03	MH0AQ5	3	1.15	9.94	8.84	87.5
K2720-05	MH0AQ9	4	1.12	9.55	8.38	86.1
K2720-06	MH0AR0	5	1.19	9.72	8.4	84.5
K2720-07	MH0AR1	6	1.13	9.94	8.74	86.4
K2720-08	MH0AR2	7	1.18	9.58	8.39	85.8
K2720-09	MH0AR3	8	1.14	9.52	8.4	86.6
K2720-10	MH0AR4	9	1.16	9.59	8.54	87.5
K2720-11	MH0AR5	10	1.18	9.94	8.58	84.5
K2720-12	MH0AR6	11	1.19	9.94	8.79	86.9
K2720-13	MH0AR7	12	1.17	9.55	8.42	86.5
K2720-14	MH0AR8	13	1.11	9.61	7.78	78.5
K2720-15	MH0AR9	14	1.18	9.93	8.00	77.9
K2720-16	MH0AS0	15	1.18	9.64	7.48	74.5
K2720-17	MH0AS1	16	1.15	9.74	8.78	88.8
K2720-18	MH0AS2	17	1.17	9.83	8.95	89.8
K2720-19	MH0AS3	18	1.19	9.94	9.15	91.0
K2720-20	MH0AS4	19	1.13	9.66	8.81	90.0
K2720-21	MH0AS4D	20	1.13	9.66	8.81	90.0
K2720-22	MH0AS4S	21	1.13	9.66	8.81	90.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$