



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:55
In Date: 05/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:34
Out Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102605

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2723-01	MH0AW2	1	1.16	9.68	8.83	90.0
K2723-02	MH0AW3	2	1.13	9.52	8.63	89.4
K2723-03	MH0AW4	3	1.18	9.55	8.71	90.0
K2723-04	MH0AW5	4	1.15	9.66	8.86	90.6
K2723-05	MH0AW6	5	1.19	9.81	8.83	88.6
K2723-06	MH0AW7	6	1.11	9.59	8.57	88.0
K2723-07	MH0AW8	7	1.14	9.65	9.22	94.9
K2723-08	MH0AW9	8	1.1	9.6	9.12	94.4
K2723-09	MH0AX0	9	1.12	9.95	9.48	94.7
K2723-10	MH0AX1	10	1.17	9.56	8.72	90.0
K2723-11	MH0AX2	11	1.17	9.67	8.89	90.8
K2723-12	MH0AX3	12	1.14	9.59	8.92	92.1
K2723-13	MH0AX4	13	1.18	9.83	9.38	94.8
K2723-14	MH0AX5	14	1.18	9.92	9.57	96.0
K2723-15	MH0AX6	15	1.18	9.9	9.61	96.7
K2723-16	MH0AX7	16	1.13	9.7	9.18	93.9
K2723-17	MH0AX8	17	1.17	9.88	9.34	93.8
K2723-18	MH0AX9	18	1.11	9.66	9.19	94.5
K2723-19	MH0AY0	19	1.15	9.55	8.78	90.8
K2723-20	MH0AY1	20	1.14	9.62	8.94	92.0
K2723-21	MH0AY1D	21	1.14	9.62	8.94	92.0
K2723-22	MH0AY1S	22	1.14	9.62	8.94	92.0

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$