



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 05/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:24
Out Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102612

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2733-01	P021-SS001-0002-01	39	1.13	9.65	7.94	79.9
K2733-02	K2733-01MS	40	1.13	9.65	7.94	79.9
K2733-03	K2733-01MSD	41	1.13	9.65	7.94	79.9
K2733-04	P021-SS001-0002-02	42	1.16	9.6	8.01	81.2
K2733-05	P021-SS001-0206-01	43	1.18	9.64	8.19	82.9
K2733-06	P021-SS001-0612-01	44	1.15	9.82	8.46	84.3
K2733-07	P021-SS001-1218-01	45	1.15	9.95	8.7	85.8
K2733-08	P021-SS001-1824-01	46	1.17	9.54	8.28	84.9
K2733-09	P021-SS002-0002-01	47	1.11	9.91	8.79	87.3
K2733-10	P021-SS002-0206-01	48	1.18	9.88	8.43	83.3
K2733-11	P021-SS002-0612-01	49	1.17	9.98	8.96	88.4
K2733-12	P021-SS002-1218-01	50	1.18	9.92	8.43	83.0
K2733-13	P021-SS002-1824-01	51	1.18	9.8	8.06	79.8
K2733-14	P021-SS003-0002-01	52	1.16	9.56	7.6	76.7
K2733-15	P021-SS003-0206-01	53	1.19	9.89	8.11	79.5
K2733-16	P021-SS003-0612-01	54	1.17	9.54	7.75	78.6
K2733-17	P021-SS003-1218-01	55	1.19	9.97	8.27	80.6
K2733-18	P021-SS003-1824-01	56	1.11	9.72	8.13	81.5
K2733-19	P021-SS004-0002-01	57	1.13	9.53	7.55	76.4
K2733-20	P021-SS004-0206-01	58	1.1	9.71	8.14	81.8
K2733-21	P021-SS004-0612-01	59	1.15	9.97	8.86	87.4
K2734-01	P021-SS004-1218-01	1	1.12	9.52	8.42	86.9
K2734-02	P021-SS004-1824-01	2	1.13	9.76	8.64	87.0
K2734-03	P026-SS012-0002-01	3	1.19	9.8	7.62	74.7
K2734-04	K2734-03MS	4	1.19	9.8	7.62	74.7
K2734-05	K2734-03MSD	5	1.19	9.8	7.62	74.7
K2734-06	P026-SS012-0002-02	6	1.13	9.84	7.72	75.7
K2734-07	P026-SS012-0206-01	7	1.16	9.9	8.68	86.0
K2734-08	P026-SS012-0612-01	8	1.18	9.5	8.51	88.1
K2734-09	P026-SS012-1218-01	9	1.13	9.96	8.72	86.0
K2734-10	P026-SS012-1824-01	10	1.18	9.78	8.27	82.4
K2734-11	P026-SS013-0002-01	11	1.13	9.6	7.42	74.3
K2734-12	P026-SS013-0206-01	12	1.17	9.67	7.77	77.6
K2734-13	P026-SS013-0612-01	13	1.16	9.61	7.74	77.9
K2734-14	P026-SS013-1218-01	14	1.16	9.98	8.49	83.1
K2734-15	P026-SS013-1824-01	15	1.1	9.89	8.47	83.8
K2734-16	P026-SS014-0002-01	16	1.12	9.76	6.21	58.9
K2734-17	P026-SS014-0206-01	17	1.19	9.63	7.21	71.3
K2734-18	P026-SS014-0612-01	18	1.18	9.82	7.83	77.0



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/9/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 05/08/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:24
Out Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102612

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2734-19	P026-SS014-1218-01	19	1.11	9.7	8.11	81.5
K2734-20	P026-SS014-1824-01	20	1.15	9.91	9.01	89.7
K2734-21	P026-SS015-0002-01	21	1.19	9.71	8.37	84.3
K2735-01	P026-SS015-0206-01	22	1.14	9.85	9.06	90.9
K2735-02	P026-SS015-0612-01	23	1.14	9.97	9.04	89.5
K2735-03	P026-SS015-1218-01	24	1.14	9.97	9.05	89.6
K2735-04	P026-SS015-1824-01	25	1.17	9.51	8.46	87.4
K2735-05	K2735-04MS	26	1.17	9.51	8.46	87.4
K2735-06	K2735-04MSD	27	1.17	9.51	8.46	87.4
K2735-07	P026-SS015-1824-02	28	1.16	9.89	8.97	89.5
K2735-08	P026-SS016-0002-01	29	1.13	9.78	8.06	80.1
K2735-09	P026-SS016-0206-01	30	1.14	9.98	8.72	85.7
K2735-10	P026-SS016-0612-01	31	1.12	9.98	8.42	82.4
K2735-11	P026-SS016-1218-01	32	1.17	9.77	8.62	86.6
K2735-12	P026-SS016-1824-01	33	1.11	9.77	8.48	85.1
K2735-13	P026-SS017-0002-01	34	1.15	9.96	5.46	48.9
K2735-14	P026-SS017-0206-01	35	1.14	9.52	7.37	74.3
K2735-15	P026-SS017-0612-01	36	1.14	9.96	7.48	71.9
K2735-16	P026-SS017-1218-01	37	1.16	9.96	8.08	78.6
K2735-17	P026-SS017-1824-01	38	1.12	9.97	7.79	75.4
K2746-01	COMP-5.8.19	60	1.12	9.6	8.89	91.6
K2749-01	TRS-1	61	1.00	2.00	2.00	100.0
K2749-03	TRS-2	62	1.00	2.00	2.00	100.0
K2749-05	FR-COMP	63	1.00	2.00	2.00	100.0
K2750-01	VB16199	64	1.19	9.88	8.16	80.2
K2750-03	72-11995	65	1.00	2.00	2.00	100.0
K2750-04	72-11995	66	1.00	2.00	2.00	100.0
K2750-05	72-11995	67	1.00	2.00	2.00	100.0
K2751-01	COMP-3	68	1.00	2.00	2.00	100.0
K2751-02	COMP-3	69	1.00	2.00	2.00	100.0
K2754-01	TP-4	70	1.12	9.55	7.86	80.0
K2754-02	TP-5	71	1.16	9.87	8.49	84.2
K2754-03	TP-3	72	1.18	9.98	8.54	83.6
K2754-04	TP-1	73	1.15	9.85	8.25	81.6
K2754-05	C-1	74	1.11	9.88	8.45	83.7
K2754-07	C-2	75	1.19	9.98	8.58	84.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$