



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:30
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:30
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102630

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2745-01	MH0H31	1	1.14	9.66	7.68	76.8
K2745-02	MH0H32	2	1.18	9.77	6.35	60.2
K2745-03	MH0H33	3	1.12	9.67	7.25	71.7
K2745-04	MH0H34	4	1.14	9.62	6.21	59.8
K2745-05	MH0H35	5	1.17	9.73	7.59	75.0
K2745-06	MH0H36	6	1.13	9.55	8.91	92.4
K2745-07	MH0H37	7	1.17	9.8	7.82	77.1
K2745-08	MH0H38	8	1.14	9.63	6.53	63.5
K2745-09	MH0H39	9	1.14	9.53	6.8	67.5
K2745-10	MH0H40	10	1.16	9.5	7.48	75.8
K2745-11	MH0H41	11	1.16	9.8	8.78	88.2
K2745-12	MH0H42	12	1.15	9.62	8.73	89.5
K2745-13	MH0H43	13	1.13	9.5	8.00	82.1
K2745-14	MH0H44	14	1.17	9.74	8.79	88.9
K2745-15	MH0H45	15	1.14	9.98	7.79	75.2
K2745-16	MH0H46	16	1.16	9.61	8.25	83.9
K2745-17	MH0H47	17	1.13	9.96	9.36	93.2
K2745-18	MH0H47D	18	1.13	9.96	9.36	93.2
K2745-19	MH0H47S	19	1.13	9.96	9.36	93.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$