



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102632

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2642-03	OK-03-050919	45	1.15	9.94	9.3	92.7
K2647-05	SU-02-050919-A	46	1.15	9.54	8.88	92.1
K2647-07	SU-02-050919-B	47	1.14	9.97	9.11	90.3
K2728-06	P-3A-OW(0-10)	1	1.11	9.91	8.73	86.6
K2756-01	P021-SS005-0002-01	2	1.12	9.96	7.93	77.0
K2756-02	P021-SS005-0206-01	3	1.14	9.65	8.29	84.0
K2756-03	P021-SS005-0612-01	4	1.16	9.53	8.73	90.4
K2756-04	P021-SS005-1218-01	5	1.16	9.77	9.05	91.6
K2756-05	P021-SS005-1824-01	6	1.18	9.56	9.09	94.4
K2756-06	P021-SS006-0002-01	7	1.14	9.69	7.5	74.4
K2756-07	P021-SS006-0206-01	8	1.16	9.77	8.42	84.3
K2756-08	P021-SS006-0612-01	9	1.18	9.72	9.12	93.0
K2756-09	P021-SS006-1218-01	10	1.12	9.79	9.08	91.8
K2756-10	P021-SS006-1824-01	11	1.18	9.71	9.13	93.2
K2756-11	P021-SS007-0002-01	12	1.12	9.78	7.89	78.2
K2756-12	P021-SS007-0206-01	13	1.18	9.98	7.92	76.6
K2756-13	P021-SS007-0612-01	14	1.13	9.57	8.05	82.0
K2756-14	P021-SS007-1218-01	15	1.14	9.66	8.24	83.3
K2756-15	P021-SS007-1824-01	16	1.17	9.96	9.11	90.3
K2756-16	P021-SS008-0002-01	17	1.15	9.92	7.58	73.3
K2756-17	P021-SS008-0206-01	18	1.18	9.74	8.16	81.5
K2756-18	K2756-17MS	19	1.18	9.74	8.16	81.5
K2756-19	K2756-17MSD	20	1.18	9.74	8.16	81.5
K2756-20	P021-SS008-0206-02	21	1.1	9.6	8.37	85.5
K2756-21	P021-SS008-0612-01	22	1.16	9.94	8.8	87.0
K2757-01	P021-SS008-1218-01	23	1.15	9.5	8.46	87.5
K2757-02	P021-SS008-1824-01	24	1.13	9.93	8.6	84.9
K2757-03	P021-SS009-0002-01	25	1.16	9.84	6.99	67.2
K2757-04	P021-SS009-0206-01	26	1.17	9.94	8.13	79.4
K2757-05	P021-SS009-0612-01	27	1.17	9.83	8.08	79.8
K2757-06	P021-SS010-0002-01	28	1.11	9.59	7.84	79.4
K2757-07	P021-SS010-0206-01	29	1.18	9.98	8.67	85.1
K2757-08	P021-SS010-0612-01	30	1.14	9.72	8.7	88.1
K2759-01	FK-0D015-03-004	31	1.14	9.52	9.22	96.4
K2759-02	FK-0D015-03-004	32	1.19	9.59	9.33	96.9
K2760-01	FK-0D015-03-005	33	1.18	9.68	8.59	87.2
K2760-02	FK-0D015-03-005	34	1.16	9.62	8.59	87.8
K2761-01	FK-SF-03-006-A	35	1.16	9.56	8.13	83.0
K2761-02	FK-SF-03-006-A	36	1.13	9.57	8.6	88.5



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:15
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:26
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102632

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2761-03	FK-SF-03-006-B	37	1.17	9.76	8.79	88.7
K2761-04	FK-SF-03-006-B	38	1.11	9.94	8.92	88.4
K2770-01	OILY-STONES	48	1.00	2.00	2.00	100.0
K2770-02	OILY-STONES	49	1.00	2.00	2.00	100.0
K2770-04	SOIL-DRUM	50	1.1	9.69	8.55	86.7
K2770-05	SOIL-DRUM	51	1.1	9.69	8.55	86.7
K2771-01	FK-PS-01-047	39	1.17	9.59	8.58	88.0
K2771-02	FK-PS-01-047	40	1.15	9.9	9.11	91.0
K2774-01	FK-PS-01-048	41	1.15	9.57	8.51	87.4
K2774-02	FK-PS-01-048	42	1.16	9.95	9.09	90.2
K2779-01	RT-5124	52	1.16	9.79	9.28	94.1
K2779-03	RBR-200050	53	1.18	9.97	8.28	80.8
K2779-05	RBR-200052	54	1.17	9.97	8.63	84.8
K2790-01	SSSI-0509-DH-4	43	1.12	9.58	8.73	90.0
K2790-02	SSSI-0509-DH-4	44	1.17	9.54	8.55	88.2
K2801-01	RT-1873	55	1.18	9.85	8.41	83.4
K2801-02	RT-1873	56	1.18	9.85	8.41	83.4
K2801-04	RBR-200035	57	1.17	9.92	9.37	93.7
K2801-05	RBR-200035	58	1.17	9.92	9.37	93.7
K2801-07	RT-3456	59	1.17	9.77	8.85	89.3
K2801-08	RT-3456	60	1.17	9.77	8.85	89.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$