



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:15
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102663

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2762-01	P021-SS001-0002-01	1	1.13	9.65	7.94	79.9
K2762-02	K2762-01MS	2	1.13	9.65	7.94	79.9
K2762-03	K2762-01MSD	3	1.13	9.65	7.94	79.9
K2762-04	P021-SS001-0002-02	4	1.16	9.6	8.01	81.2
K2762-05	P021-SS001-0206-01	5	1.18	9.64	8.19	82.9
K2762-06	P021-SS001-0612-01	6	1.15	9.82	8.46	84.3
K2762-07	P021-SS001-1218-01	7	1.15	9.95	8.7	85.8
K2762-08	P021-SS001-1824-01	8	1.17	9.54	8.28	84.9
K2762-09	P021-SS002-0002-01	9	1.11	9.91	8.79	87.3
K2762-10	P021-SS002-0206-01	10	1.18	9.88	8.43	83.3
K2762-11	P021-SS002-0612-01	11	1.17	9.98	8.96	88.4
K2762-12	P021-SS002-1218-01	12	1.18	9.92	8.43	83.0
K2762-13	P021-SS002-1824-01	13	1.18	9.8	8.06	79.8
K2762-14	P021-SS003-0002-01	14	1.16	9.56	7.6	76.7
K2762-15	P021-SS003-0206-01	15	1.19	9.89	8.11	79.5
K2762-16	P021-SS003-0612-01	16	1.17	9.54	7.75	78.6
K2762-17	P021-SS003-1218-01	17	1.19	9.97	8.27	80.6
K2762-18	P021-SS003-1824-01	18	1.11	9.72	8.13	81.5
K2762-19	P021-SS004-0002-01	19	1.13	9.53	7.55	76.4
K2762-20	P021-SS004-0206-01	20	1.11	9.71	8.14	81.7
K2762-21	P021-SS004-0612-01	21	1.15	9.97	8.86	87.4
K2763-01	P021-SS004-1218-01	22	1.12	9.52	8.42	86.9
K2763-02	P021-SS004-1824-01	23	1.13	9.76	8.64	87.0
K2763-03	P021-SS005-0002-01	24	1.12	9.96	7.93	77.0
K2763-04	P021-SS005-0206-01	25	1.14	9.65	8.29	84.0
K2763-05	P021-SS005-0612-01	26	1.16	9.53	8.73	90.4
K2763-06	P021-SS005-1218-01	27	1.16	9.77	9.05	91.6
K2763-07	P021-SS005-1824-01	28	1.18	9.56	9.09	94.4
K2763-08	P021-SS006-0002-01	29	1.14	9.69	7.5	74.4
K2763-09	P021-SS006-0206-01	30	1.16	9.77	8.42	84.3
K2763-10	P021-SS006-0612-01	31	1.18	9.72	9.12	93.0
K2763-11	P021-SS006-1218-01	32	1.12	9.79	9.08	91.8
K2763-12	P021-SS006-1824-01	33	1.18	9.71	9.13	93.2
K2763-13	P021-SS007-0002-01	34	1.12	9.78	7.89	78.2
K2763-14	P021-SS007-0206-01	35	1.18	9.98	7.92	76.6
K2763-15	P021-SS007-0612-01	36	1.13	9.57	8.05	82.0
K2763-16	P021-SS007-1218-01	37	1.14	9.66	8.24	83.3
K2763-17	P021-SS007-1824-01	38	1.17	9.96	9.11	90.3
K2763-18	P021-SS008-0002-01	39	1.15	9.92	7.58	73.3



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:15
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102663

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2763-19	P021-SS008-0206-01	40	1.18	9.74	8.16	81.5
K2763-20	K2763-19MS	41	1.18	9.74	8.16	81.5
K2763-21	K2763-19MSD	42	1.18	9.74	8.16	81.5
K2764-01	P021-SS008-0206-02	43	1.1	9.6	8.37	85.5
K2764-02	P021-SS008-0612-01	44	1.16	9.94	8.8	87.0
K2764-03	P021-SS008-1218-01	45	1.15	9.5	8.46	87.5
K2764-04	P021-SS008-1824-01	46	1.13	9.93	8.6	84.9
K2764-05	P021-SS009-0002-01	47	1.16	9.84	6.99	67.2
K2764-06	P021-SS009-0206-01	48	1.17	9.94	8.13	79.4
K2764-07	P021-SS009-0612-01	49	1.17	9.83	8.08	79.8
K2764-08	P021-SS010-0002-01	50	1.11	9.59	7.84	79.4
K2764-09	P021-SS010-0206-01	51	1.18	9.98	8.67	85.1
K2764-10	P021-SS010-0612-01	52	1.14	9.72	8.7	88.1
K2764-11	P026-SS001-0002-01	53	1.15	9.96	8.65	85.1
K2764-12	P026-SS001-0206-01	54	1.16	9.98	9.07	89.7
K2764-13	P026-SS001-0612-01	55	1.17	9.68	9.02	92.2
K2764-14	P026-SS001-1218-01	56	1.16	9.53	8.64	89.4
K2764-15	P026-SS001-1824-01	57	1.16	9.76	8.75	88.3
K2764-16	P026-SS002-0002-01	58	1.16	9.65	7.25	71.7
K2764-17	P026-SS002-0206-01	59	1.16	9.61	7.7	77.4
K2764-18	P026-SS002-0612-01	60	1.14	9.78	8.04	79.9
K2764-19	P026-SS002-1218-01	61	1.15	9.59	7.92	80.2
K2764-20	K2764-19MS	62	1.15	9.59	7.92	80.2
K2764-21	K2764-19MSD	63	1.15	9.59	7.92	80.2
K2765-01	P026-SS002-1218-02	64	1.14	9.7	8.02	80.4
K2765-02	P026-SS002-1824-01	65	1.15	9.93	8.54	84.2
K2765-03	P026-SS003-0002-01	66	1.15	9.73	7.14	69.8
K2765-04	P026-SS003-0206-01	67	1.14	9.82	7.83	77.1
K2765-05	P026-SS003-0612-01	68	1.17	9.51	7.51	76.0
K2765-06	P026-SS003-1218-01	69	1.15	9.56	7.56	76.2
K2765-07	P026-SS003-1824-01	70	1.14	9.92	7.92	77.2
K2765-08	P026-SS004-0002-01	71	1.15	9.68	8.22	82.9
K2765-09	P026-SS004-0206-01	72	1.14	9.67	8.82	90.0
K2765-10	P026-SS004-0612-01	73	1.14	9.69	9.05	92.5
K2765-11	P026-SS004-1218-01	74	1.14	9.81	8.82	88.6
K2765-12	P026-SS004-1824-01	75	1.15	9.87	8.78	87.5
K2765-13	P026-SS005-0002-01	76	1.14	9.96	8.14	79.4
K2765-14	P026-SS005-0206-01	77	1.15	9.97	8.63	84.8
K2765-15	P026-SS005-0612-01	78	1.15	9.66	8.16	82.4



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:15
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102663

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2765-16	P026-SS005-1218-01	79	1.15	9.62	8.38	85.4
K2765-17	P026-SS005-1824-01	80	1.16	9.52	8.14	83.5
K2765-18	P026-SS006-0002-01	81	1.16	9.94	7.83	76.0
K2765-19	K2765-18MS	82	1.16	9.94	7.83	76.0
K2765-20	K2765-18MSD	83	1.16	9.94	7.83	76.0
K2765-21	P026-SS006-0002-02	84	1.15	9.92	8.23	80.7
K2766-01	P026-SS006-0206-01	85	1.14	9.93	9.14	91.0
K2766-02	P026-SS006-0612-01	86	1.16	9.97	9.4	93.5
K2766-03	P026-SS006-1218-01	87	1.14	9.94	9.06	90.0
K2766-04	P026-SS006-1824-01	88	1.14	9.57	8.55	87.9
K2766-05	P026-SS007-0002-01	89	1.13	9.98	8.06	78.3
K2766-06	P026-SS007-0206-01	90	1.12	9.97	8.09	78.8
K2766-07	P026-SS007-0612-01	91	1.16	9.74	8.11	81.0
K2766-08	P026-SS007-1218-01	92	1.15	9.51	8.00	81.9
K2766-09	P026-SS007-1824-01	93	1.14	9.93	8.38	82.4
K2766-10	P026-SS008-0002-01	94	1.15	9.85	5.44	49.3
K2766-11	P026-SS008-0206-01	95	1.14	9.77	4.81	42.5
K2766-12	K2766-11MS	96	1.14	9.77	4.81	42.5
K2766-13	K2766-11MSD	97	1.14	9.77	4.81	42.5
K2766-14	P026-SS008-0206-02	98	1.14	9.92	5.09	45.0
K2766-15	P026-SS008-0612-01	99	1.13	9.63	7.97	80.5
K2766-16	P026-SS008-1218-01	100	1.15	9.95	8.76	86.5
K2766-17	P026-SS008-1824-01	101	1.15	9.84	8.58	85.5
K2766-18	P026-SS009-0002-01	102	1.15	9.68	8.45	85.6
K2766-19	P026-SS009-0206-01	103	1.14	9.98	8.98	88.7
K2766-20	P026-SS009-0612-01	104	1.15	9.57	8.69	89.5
K2766-21	P026-SS009-1218-01	105	1.15	9.97	8.94	88.3
K2767-01	P026-SS009-1824-01	106	1.15	9.76	8.56	86.1
K2767-02	P026-SS010-0002-01	107	1.16	9.86	7.45	72.3
K2767-03	P026-SS010-0206-01	108	1.14	9.54	7.89	80.4
K2767-04	P026-SS010-0612-01	109	1.15	9.73	8.08	80.8
K2767-05	P026-SS010-1218-01	110	1.15	9.93	8.34	81.9
K2767-06	P026-SS010-1824-01	111	1.16	9.89	8.42	83.2
K2767-07	P026-SS011-0002-01	112	1.15	9.94	4.93	43.0
K2767-08	P026-SS011-0206-01	113	1.15	9.93	4.28	35.6
K2767-09	P026-SS011-0612-01	114	1.15	9.56	6.72	66.2
K2767-10	P026-SS011-1218-01	115	1.15	9.78	7.25	70.7
K2767-11	P026-SS011-1824-01	116	1.15	9.63	7.76	77.9
K2767-12	P026-SS012-0002-01	117	1.19	9.8	7.62	74.7



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/10/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 17:30
In Date: 05/09/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 102
Time OUT: 08:15
Out Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102663

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2767-13	K2767-12MS	118	1.19	9.8	7.62	74.7
K2767-14	K2767-12MSD	119	1.19	9.8	7.62	74.7
K2767-15	P026-SS012-0002-02	120	1.13	9.84	7.72	75.7
K2767-16	P026-SS012-0206-01	121	1.16	9.9	8.68	86.0
K2767-17	P026-SS012-0612-01	122	1.18	9.5	8.51	88.1
K2767-18	P026-SS012-1218-01	123	1.13	9.96	8.72	86.0
K2767-19	P026-SS012-1824-01	124	1.18	9.78	8.27	82.4
K2767-20	P026-SS013-0002-01	125	1.13	9.6	7.42	74.3
K2767-21	P026-SS013-0206-01	126	1.17	9.67	7.77	77.6
K2768-01	P026-SS013-0612-01	127	1.16	9.61	7.74	77.9
K2768-02	P026-SS013-1218-01	128	1.16	9.98	8.49	83.1
K2768-03	P026-SS013-1824-01	129	1.1	9.89	8.47	83.8
K2768-04	P026-SS014-0002-01	130	1.12	9.76	6.21	58.9
K2768-05	P026-SS014-0206-01	131	1.19	9.63	7.21	71.3
K2768-06	P026-SS014-0612-01	132	1.18	9.82	7.83	77.0
K2768-07	P026-SS014-1218-01	133	1.11	9.7	8.11	81.5
K2768-08	P026-SS014-1824-01	134	1.15	9.91	9.01	89.7
K2768-09	P026-SS015-0002-01	135	1.19	9.71	8.37	84.3
K2768-10	P026-SS015-0206-01	136	1.14	9.85	9.06	90.9
K2768-11	P026-SS015-0612-01	137	1.14	9.97	9.04	89.5
K2768-12	P026-SS015-1218-01	138	1.14	9.97	9.05	89.6
K2768-13	P026-SS015-1824-01	139	1.17	9.51	8.46	87.4
K2768-14	K2768-13MS	140	1.17	9.51	8.46	87.4
K2768-15	K2768-13MSD	141	1.17	9.51	8.46	87.4
K2768-16	P026-SS015-1824-02	142	1.16	9.89	8.97	89.5
K2768-17	P026-SS016-0002-01	143	1.13	9.78	8.06	80.1
K2768-18	P026-SS016-0206-01	144	1.14	9.98	8.72	85.7
K2768-19	P026-SS016-0612-01	145	1.12	9.98	8.42	82.4
K2768-20	P026-SS016-1218-01	146	1.17	9.77	8.62	86.6
K2768-21	P026-SS016-1824-01	147	1.11	9.77	8.48	85.1
K2769-01	P026-SS017-0002-01	148	1.15	9.96	5.46	48.9
K2769-02	P026-SS017-0206-01	149	1.14	9.52	7.37	74.3
K2769-03	P026-SS017-0612-01	150	1.14	9.96	7.48	71.9
K2769-04	P026-SS017-1218-01	151	1.16	9.96	8.08	78.6
K2769-05	P026-SS017-1824-01	152	1.12	9.97	7.79	75.4

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$