



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/13/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:00
In Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:35
Out Date: 05/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102666

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2737-01	DB848	1	1.15	9.63	9.62	99.9
K2737-02	DB850	2	1.19	9.75	9.33	95.1
K2737-03	DB850MS	3	1.19	9.75	9.33	95.1
K2737-04	DB850MSD	4	1.19	9.75	9.33	95.1
K2737-05	DB851	5	1.19	9.76	8.79	88.7
K2737-06	DB856	6	1.14	9.62	8.75	89.7
K2737-07	DB857	7	1.17	9.66	8.81	90.0
K2737-08	VHBLK01	8	1.00	2.00	2.00	100.0
K2737-09	DB847	9	1.00	2.00	2.00	100.0
K2737-10	DB858	10	1.11	9.52	9.51	99.9
K2737-11	DB859	11	1.18	9.81	9.15	92.4
K2737-12	DB860	12	1.17	9.75	9.37	95.6
K2737-13	DB861	13	1.14	9.65	8.56	87.2
K2737-14	DB862	14	1.14	9.96	9.28	92.3
K2737-15	DB863	15	1.13	9.59	9.04	93.5
K2737-16	DB864	16	1.15	9.93	9.46	94.6
K2737-17	DB865	17	1.18	9.58	9.13	94.6
K2737-18	DB866	18	1.15	9.95	9.54	95.3
K2737-19	DB867	19	1.17	9.88	9.07	90.7
K2737-20	DB868	20	1.12	9.62	9.21	95.2
K2737-21	DB869	21	1.14	9.96	9.44	94.1
K2737-22	DB870	22	1.17	9.9	9.28	92.9
K2737-23	DB871	23	1.11	9.69	8.87	90.4
K2737-24	DB872	24	1.19	9.64	8.57	87.3

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$