



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/13/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 15:35
In Date: 05/10/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 05/11/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102668

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2794-01	MGAKL2	1	1.15	9.58	7.58	76.3
K2794-02	MGAKL3	2	1.1	9.91	8.34	82.2
K2794-03	MGAKL3D	3	1.1	9.91	8.34	82.2
K2794-04	MGAKL3S	4	1.1	9.91	8.34	82.2
K2794-05	MGAKL4	5	1.14	9.67	8.09	81.5
K2794-06	MGAKL5	6	1.13	9.85	9.06	90.9
K2794-07	MGAKL6	7	1.13	9.53	6.34	62.0
K2794-08	MGAKL7	8	1.18	9.95	9.11	90.4
K2794-09	MGAKL8	9	1.11	9.59	9.3	96.6
K2794-10	MGAKL9	10	1.13	9.74	9.03	91.8
K2794-11	MGAKM0	11	1.18	9.97	9.2	91.2
K2794-12	MGAKM1	12	1.13	9.93	7.59	73.4
K2794-13	MGAKM2	13	1.16	9.87	8.12	79.9
K2794-14	MGAKM3	14	1.16	9.55	7.11	70.9
K2794-15	MGAKM4	15	1.18	9.91	7.38	71.0
K2794-16	MGAKM5	16	1.11	9.88	7.99	78.4
K2794-17	MGAKM6	17	1.13	9.66	7.62	76.1
K2794-18	MGAKM7	18	1.11	9.63	7.46	74.5
K2794-19	MGAKM8	19	1.18	9.56	7.26	72.6
K2794-20	MGAKM9	20	1.19	9.88	7.48	72.4
K2794-21	MGAKN0	21	1.14	9.52	7.16	71.8

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$