



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/14/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 17:00
In Date: 05/13/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:00
Out Date: 05/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102704

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2831-01	MA5108	1	1.00	2.00	2.00	100.0
K2831-02	MA5112	2	1.00	2.00	2.00	100.0
K2831-03	MA5115	3	1.00	2.00	2.00	100.0
K2831-04	MA5172	4	1.14	9.55	9.09	94.5
K2831-05	MA5173	5	1.16	9.85	8.86	88.6
K2831-06	MA5174	6	1.12	9.59	8.81	90.8
K2831-07	MA5175	7	1.14	9.93	8.69	85.9
K2831-08	MA5176	8	1.14	9.6	8.5	87.0
K2831-09	MA5177	9	1.18	9.54	9.02	93.8
K2831-10	MA5178	10	1.15	9.94	9.24	92.0
K2831-11	MA5179	11	1.17	9.97	9.18	91.0
K2831-12	MA5180	12	1.1	9.65	8.95	91.8
K2831-13	MA5181	13	1.13	9.79	9.33	94.7
K2831-14	MA5182	14	1.11	9.85	9.31	93.8
K2831-15	MA5183	15	1.18	9.98	9.37	93.1
K2831-16	MA5184	16	1.17	9.7	9.22	94.4
K2831-17	MA5185	17	1.15	9.73	8.68	87.8
K2831-18	MA5186	18	1.14	9.96	8.77	86.5
K2831-19	MA5186D	19	1.14	9.96	8.77	86.5
K2831-20	MA5186S	20	1.14	9.96	8.77	86.5
K2831-21	MA5187	21	1.11	9.95	8.93	88.5
K2831-22	MA5188	22	1.16	9.76	8.67	87.3
K2831-23	MA5189	23	1.15	9.95	9.45	94.3
K2831-24	MA5190	24	1.19	9.86	9.00	90.1
K2831-25	MA5191	25	1.19	9.86	9.22	92.6

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$