



PERCENT SOLID

Supervisor: apatel
Analyst: jignesh
Date: 5/15/2019

OVENTEMP IN Celsius(°C): 107
Time IN: 16:00
In Date: 05/14/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:45
Out Date: 05/15/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102738

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2832-01	MA5196	1	1.15	9.96	9.25	91.9
K2832-02	MA5197	2	1.14	9.82	8.28	82.3
K2832-03	MA5198	3	1.19	9.58	9.00	93.1
K2832-04	MA5199	4	1.13	9.73	9.14	93.1
K2832-05	MA51A0	5	1.18	9.75	8.49	85.3
K2832-06	MA51A1	6	1.14	9.92	9.02	89.7
K2832-08	MA51A3	7	1.12	9.57	9.04	93.7
K2832-09	MA51A4	8	1.16	9.97	8.92	88.1
K2832-10	MA51A5	9	1.15	9.91	8.98	89.4
K2832-11	MA51A6	10	1.11	9.76	8.59	86.5
K2832-12	MA5192	11	1.11	9.7	8.71	88.5
K2832-13	MA5193	12	1.16	9.76	8.98	90.9
K2832-14	MA51A8	13	1.17	9.88	8.64	85.8
K2832-15	MA51A9	14	1.12	9.83	8.24	81.7
K2832-16	MA51B0	15	1.18	9.85	8.69	86.6
K2832-17	MA51B0D	16	1.18	9.85	8.69	86.6
K2832-18	MA51B0S	17	1.18	9.85	8.69	86.6
K2832-19	MA51B1	18	1.13	9.68	8.53	86.5

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$