

Test results

Aquakem 7.2AQ1

Page: 1

CHEMTECH
284 Sheffield Street,
Mountainside, NJ 07092
Reviewed by : AS

Instrument ID : Konelab

5/15/2019 13:36

Test: ~~Total~~ ^{Reactive} CN

Sample Id	Result ^{ug/d}	Dil. 1 +	Response	Errors
ICV1	90.096	0.0	0.093	
ICB1	0.200	0.0	0.004	
CCV1	242.437	0.0	0.244	
CCB1	-0.187	0.0	0.004	
PB119725BL	-0.387	0.0	0.004	
K2836-03	-0.404	0.0	0.004	
K2836-03DUP	-0.555	0.0	0.004	
K2837-02	-0.276	0.0	0.004	
PB119756BL	-1.007	0.0	0.003	
K2852-02	-0.234	0.0	0.004	
K2852-02DUP	-0.166	0.0	0.004	
PB119755BL	-0.486	0.0	0.004	
K2851-03	-0.429	0.0	0.004	
K2851-03DUP	-0.572	0.0	0.004	
CCV2	256.291	0.0	0.258	
CCB2	-0.155	0.0	0.004	

N 16
Mean 36.510
SD 86.1211
CV% 235.88

Aquakem v. 7.2AQ1

Results from time period:

Wed May 15 12:57:25 2019

Wed May 15 13:33:05 2019

AS

Sample Id	Sam/Ctr/	Test short nam	Test type	Result	Result unit	Result date and time	Stat
0.OPPBCN	A	Reactive CN	P	-0.6691	µg/l	5/15/2019 12:57:25	
5.OPPBCN	A	Reactive CN	P	4.4965	µg/l	5/15/2019 12:57:26	
10PPBCN	A	Reactive CN	P	9.902	µg/l	5/15/2019 12:57:27	
50PPBCN	A	Reactive CN	P	49.5709	µg/l	5/15/2019 12:57:28	
100PPBCN	A	Reactive CN	P	99.7831	µg/l	5/15/2019 12:57:29	
250PPBCN	A	Reactive CN	P	253.6466	µg/l	5/15/2019 12:57:30	
500PPBCN	A	Reactive CN	P	498.27	µg/l	5/15/2019 12:57:31	
ICV1	S	Reactive CN	P	90.0962	µg/l	5/15/2019 13:28:14	
ICB1	S	Reactive CN	P	0.2004	µg/l	5/15/2019 13:28:15	
CCV1	S	Reactive CN	P	242.4372	µg/l	5/15/2019 13:28:16	
CCB1	S	Reactive CN	P	-0.187	µg/l	5/15/2019 13:28:17	
PB119725BL	S	Reactive CN	P	-0.3871	µg/l	5/15/2019 13:28:18	
K2836-03	S	Reactive CN	P	-0.4043	µg/l	5/15/2019 13:28:19	
K2836-03DUP	S	Reactive CN	P	-0.5548	µg/l	5/15/2019 13:28:20	
K2837-02	S	Reactive CN	P	-0.2757	µg/l	5/15/2019 13:28:21	
PB119756BL	S	Reactive CN	P	-1.0074	µg/l	5/15/2019 13:28:22	
K2852-02	S	Reactive CN	P	-0.2338	µg/l	5/15/2019 13:28:23	
K2852-02DUP	S	Reactive CN	P	-0.1662	µg/l	5/15/2019 13:28:24	
PB119755BL	S	Reactive CN	P	-0.4855	µg/l	5/15/2019 13:33:01	
K2851-03	S	Reactive CN	P	-0.4292	µg/l	5/15/2019 13:33:02	
K2851-03DUP	S	Reactive CN	P	-0.5723	µg/l	5/15/2019 13:33:03	
CCV2	S	Reactive CN	P	256.2907	µg/l	5/15/2019 13:33:04	
CCB2	S	Reactive CN	P	-0.155	µg/l	5/15/2019 13:33:05	

CHEMTECH
 284 Sheffield Street,
 Mountainside, NJ 07092
 Reviewed by : AS

Instrument ID : Konelab

5/15/2019 13:04


Test Total CN

Accepted 5/15/2019 13:04

Factor 1009
 Bias 0.004

Coeff. of det. 0.999916

Errors



	Calibrator	Response	Calc. con.	Conc.	Errors
1	0.0PPBCN	0.003	-0.6691	0.0000	
2	5.0PPBCN	0.009	4.4965	5.0000	
3	10PPBCN	0.014	9.9020	10.0000	
4	50PPBCN	0.053	49.5709	50.0000	
5	100PPBCN	0.103	99.7831	100.0000	
6	250PPBCN	0.256	253.6466	250.0000	
7	500PPBCN	0.498	498.2700	500.0000	