



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/20/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 16:45
In Date: 05/17/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 08:15
Out Date: 05/18/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102828

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2925-01	A4284	1	1.15	9.95	9.4	93.8
K2925-02	A4285	2	1.16	9.58	8.68	89.3
K2925-03	A4286	3	1.19	9.69	8.43	85.2
K2925-04	A4287	4	1.11	9.55	7.36	74.1
K2925-05	A4288	5	1.18	9.81	7.74	76.0
K2925-06	A4288MS	6	1.18	9.81	7.74	76.0
K2925-07	A4288MSD	7	1.18	9.81	7.74	76.0
K2925-08	A4289	8	1.15	9.54	8.34	85.7
K2925-09	A4290	9	1.17	9.92	9.45	94.6
K2925-10	A4291	10	1.18	9.98	9.38	93.2
K2925-11	A4296	11	1.14	9.94	7.29	69.9
K2925-12	A4297	12	1.17	9.9	7.15	68.5
K2925-13	A4298	13	1.16	9.93	7.49	72.2
K2925-14	A4299	14	1.15	9.95	8.34	81.7
K2925-15	A42A0	15	1.13	9.93	7.51	72.5
K2925-16	A42A1	16	1.18	9.5	7.08	70.9
K2925-17	A42A2	17	1.17	9.89	7.88	76.9
K2925-18	A42A3	18	1.15	9.84	8.06	79.5
K2925-20	A42E9	19	1.00	2.00	2.00	100.0
K2925-21	A42F5	20	1.19	9.7	7.24	71.1

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$