



PERCENT SOLID

Supervisor: apatel
Analyst: JIGNESH
Date: 5/20/2019

OVENTEMP IN Celsius(°C): 108
Time IN: 14:00
In Date: 05/19/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
OvenID: M OVEN-1

OVENTEMP OUT Celsius(°C): 103
Time OUT: 07:20
Out Date: 05/20/2019
Weight Check 1.0g: 1.00
Weight Check 10g: 10.00
BalanceID: M SC-1
Thermometer ID: %SOLIDS-OVEN

QC:LB102844

Lab ID	Client Sample ID	Dish#	Dish Wt(g) (A)	Dish + Sample Wt(g) (B)	Dish + Dry Sample Wt(g) (C)	% Solid
K2927-01	A4292	1	1.18	9.84	8.32	82.4
K2927-02	A4293	2	1.16	9.92	8.52	84.0
K2927-03	A4294	3	1.19	9.93	8.84	87.5
K2927-04	A4295	4	1.16	9.72	8.84	89.7
K2927-05	A42A4	5	1.15	9.96	8.62	84.8
K2927-06	A42A4MS	6	1.15	9.96	8.62	84.8
K2927-07	A42A4MSD	7	1.15	9.96	8.62	84.8
K2927-08	A42A5	8	1.17	9.72	8.08	80.8
K2927-09	A42A6	9	1.14	9.85	9.13	91.7
K2927-10	A42A7	10	1.18	9.87	8.99	89.9
K2927-11	A42A8	11	1.16	9.96	8.4	82.3
K2927-12	A42A9	12	1.1	9.73	8.21	82.4
K2927-13	A42B0	13	1.12	9.67	8.5	86.3
K2927-14	A42B1	14	1.18	9.97	8.79	86.6
K2927-15	A42B2	15	1.14	9.52	7.54	76.4
K2927-16	A42B3	16	1.18	9.64	7.98	80.4
K2927-17	A42B4	17	1.12	9.59	7.94	80.5
K2927-18	A42B5	18	1.16	9.68	8.63	87.7
K2927-19	A42C0	19	1.14	9.94	5.56	50.2
K2927-21	A42F0	20	1.00	2.00	2.00	100.0
K2927-22	A42F6	21	1.15	9.76	8.23	82.2

$$\% \text{ Solid} = \frac{(C-A) * 100}{(B-A)}$$